JOUAN

**CENTRIFUGEURS
C3.12 / CR3.12
CENTRIFUGES**

**MANUEL D'UTILISATION
USER'S MANUAL**

25990017-h

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GUARANTEE TERMS

JOUAN guarantees that this unit is free from defects in materials and workmanship when it leaves the factory, and undertakes to replace or repair the unit if it proves defective in normal use or during servicing for a period of **ONE YEAR** from the delivery.

Our liability under this guarantee is limited to repairing the defective unit or any part of the unit providing it is sent, carriage paid, to an authorised service centre or the SAINT-HERBLAIN office.

This guarantee is invalidated if the unit is incorrectly used, poorly serviced or neglected, mis-used or accidentally damaged.

There is no explicit guarantee other than as stated above.

FOR FURTHER INFORMATION, ASSISTANCE OR SERVICING

Please contact :

USA : WINCHESTER, VIRGINIA Tel. (540) 869 8623 - Fax. (540) 869 8626

ENGLAND : ILKESTON, DERBYS Tel. (0115) 944 7989 - Fax. (0115) 944 7080

GERMANY : UNTERHACHING Tel. (089) 611 4038 - Fax. (089) 611 3087

OTHER COUNTRIES : Your local representative

or our CUSTOMER SERVICE CENTRE :
SOCIETE JOUAN
Service Assistance Technique
10, rue Duguay Trouin
44800 SAINT HERBLAIN
Tél. +33 (0) 2 28 03 20 20
Fax. +33 (0) 2 28 03 20 03

PACKING LIST

QUANTITY	ITEM NAME
1	Centrifuge
1	User's manual
1	Unpacking instructions
1	Imbalance protection instructions
1	Rotor removal tool
1	Mains cord
1	Box wrench
1	12 and 14 mm wrench
1	Allen key
1	Screwdriver
1	Screwdriver 2/35
1	Fuse
2	Carbon brushes

CONTENTS

CHAPTER 1 - USE AND FUNCTION.....	1-1
1.1 - Description	1-1
1.2 - Relative Centrifugal Force.....	1-2
CHAPTER 2 - INSTALLATION PROCEDURE.....	2-1
2.1 - Lifting and transport.....	2-1
2.2 - Unpacking	2-1
2.3 - Environmental conditions.....	2-1
2.4 - Location	2-1
2.5 - Mains supply	2-2
2.6 - Lid opening	2-2
2.7 - Manual lid opening.....	2-2
2.8 - Accessory installation.....	2-3
2.9 - Balancing the load	2-3
CHAPTER 3 - SPECIFICATIONS.....	3-1
3.1 - Dimensions and weight	3-1
3.2 - Electrical specifications	3-1
3.3 - Performance and accessories.....	3-1
CHAPTER 4 - OPERATING PRINCIPLES.....	4-1
4.1 - Drive system	4-1
4.2 - Safety interlock system	4-1
4.3 - Imbalance detection system.....	4-1
4.4 - Refrigeration system.....	4-1
CHAPTER 5 - INSTRUCTIONS FOR USE.....	5-1
5.1 - Controls and indicators	5-1
5.2 - Preliminary operations	5-2
5.3 - Centrifugation cycle	5-2
5.4 - End of the cycle	5-2
CHAPTER 6 - HAZARDS, PRECAUTIONS AND LIMITATIONS OF USE.....	6-1
6.1 - Speed control.....	6-1
6.2 - Cautions	6-1
6.3 - IEC 1010-2-020 standard	6-1
6.4 - Operational limitations.....	6-1
6.5 - Samples	6-2
6.6 - Hazards.....	6-2
6.7 - Improper use	6-2
6.8 - Biological risks	6-2
CHAPTER 7 - SERVICING AND PREVENTATIVE MAINTENANCE.....	7-1
7.1 - Cleaning	7-1
7.2 - Maintenance.....	7-1
7.3 - Rotor removal	7-2
7.4 - Fuse replacement.....	7-2
7.5 - Spare part list	7-2

CHAPTER 1

USE AND FUNCTION

1.1 - DESCRIPTION

C3.12 / CR3.12 centrifuges are designed for use in the laboratory for the separation of substances with varying density, by centrifugal force.

The CR3.12 is a refrigerated version of the C3.12 ventilated centrifuge. The refrigeration is used to maintain sample at, or below, ambient temperature.


Figure 1.1 - GENERAL VIEW

1.2 - RELATIVE CENTRIFUGAL FORCE

The C3.12 and CR3.12 are automated centrifuges designed for general purpose clinical and research laboratory use that requires the application of relative centrifugal force (R.C.F.).

Suspended solids are separated by the Relative Centrifugal Force (RCF) applied to the sample. The effective force increases with the square of the speed of rotation and the distance from the axis of rotation.

The following formulae permit the calculation of primary parameters and of transformations relating to changes in primary parameters.

In fact, the control system of the instrument carries out and displays the results of all calculations between speed, radius and RCF.

Note : The value introduced for the radius can be adjusted to allow for position within the tube such as at a boundary. Maximum radii are quoted in the specifications tables (chapter 3).

CENTRIFUGATION FORMULAE

Legend : R = radius (in mm) N = speed (in r.p.m.) ÷ 1000 RCF = gravitational acceleration 'g'
M+ = add to memory MR = memory recall

Note : To calculate actual results, press the keys on a pocket calculator in the order shown.

Primary calculations

$$\text{RCF (x g)} = 1.118 R N^2$$

$$\text{Speed (r.p.m.)} = 946 \sqrt{\frac{\text{RCF}}{R}}$$

$$\text{Radius (mm)} = \frac{\text{RCF}}{1.118 R N^2}$$

Key sequence

$$N \times = \times 1.118 \times R =$$

$$\text{RCF} \div R = \sqrt{\times 946} =$$

$$N \times = \times 1.118 = \text{M+ RCF} \div \text{MR} =$$

Transformations

To determine actual 'g' achieved at a different speed :

$$\text{RCF2} = \text{RCF1} \left(\frac{N2}{N1} \right)^2$$

$$N2 \div N1 = \times \times \text{RCF1} =$$

To determine actual speed required to achieve a different 'g' at the same radius :

$$N2 = N1 \sqrt{\frac{\text{RCF2}}{\text{RCF1}}}$$

$$\text{RCF2} \div \text{RCF1} = \sqrt{\times N1} =$$

To determine actual speed required to achieve the same 'g' at a different radius :

$$N2 = N1 \sqrt{\frac{R1}{R2}}$$

$$R1 \div R2 = \sqrt{\times N1} =$$

CHAPTER 2

INSTALLATION PROCEDURE

2.1 - LIFTING AND TRANSPORT

Because of the weight of the machine, all lifting and transporting must be carried out using proper handling equipment (eg : fork lift trolley) that complies with current regulations, and by people having undergone the necessary training.

The machine must be supported from underneath. If it has to be transported without its pallet, for example on a staircase, professional handling assistance is required.

2.2 - UNPACKING

Carefully remove the unit, the accessories and the documents from the packing.

A minimum of two people are required to handle the machine.

Check all items are present using the delivery note.

Before connecting the machine to the mains, remove the foam pad locking the imbalance detector according to the specific instruction sheet included in the unit.

2.3 - ENVIRONMENTAL CONDITIONS

This instrument is designed to operate safely under the following conditions :

- Indoor use.
- Temperature : 5°C to 40°C.
- Maximum relative humidity of 80% for temperatures up to 22°C.
- Maximum altitude : 2000 m.

Maximum performance is assured across the following temperature range : 15°C to 25°C.

2.4 - LOCATION

Locate the centrifuge on a bench able to support its weight and vibrations, in a clean, non corrosive environment.

Allow 15 cm clearance all around the centrifuge so that correct ventilation is obtained (see also IEC 1010-2-020 recommendations, chapter 6.3).

2.5 - MAINS SUPPLY

C/CR3.12	230 V \pm 10% - 50 Hz	6 A	simple phase + ground
	120 V \pm 10% - 60 Hz	10 A	simple phase + ground

Remember that in order to respect the electrical safety standards related to protection of operators against indirect contact, the supply of power to the instrument must be via a power socket fitted with a protection device ensuring automatic cut-off in the case of an insulation fault.

A supply fitted with a circuit breaker of the correct rating complies with this requirement.

2.6 - LID OPENING

Ensure the centrifuge has been switched ON. Press the button on the right hand side of the centrifuge and hold it in while pressing the LID key : the lid is automatically unlocked and opens.

Wait until speed display shows "0 rpm" before pressing the LID key.

2.7 - MANUAL LID OPENING

In the case of a mains power cut, opening of the lid is prevented by the lid locking safety device. It is recommended to wait for the mains to be switched back on so that this safety device enables the lid to be unlocked.

During rotation, the slightest leak from a sample is enough to create an aerosol. Rotation of the rotor can only be detected when the machine is switched on. Opening the lid manually if the rotor is still rotating, would disperse the aerosol in the environment.

Rotating parts are also a risk as these could come in contact with the user or be ejected. There is particularly high risk of injury if :

- The user tries to manually stop the rotor,
- Any object falls inside the centrifugation chamber.

MANUAL LID OPENING PROCEDURE

Manual lid unlocking must only be done by someone informed of the danger and of the necessary precautions to be undertaken.

If the machine has stopped because of a mains power cut, the rotor could still be rotating. Wait 10 minutes before opening the lid. In spite of the absence of noise, the rotor could still be rotating when you need to open the lid manually.

Always set the mains switch to the OFF position, even in the case of a mains power cut.

- Insert the unlocking tool into the hole at the top on the left hand side of the instrument by piercing the protection sticker.
- Insert the tool and push to unlock the lid.

The protection sticker must be replaced after it has been pierced.

2.8 - ACCESSORY INSTALLATION

- Remove the shaft nut (clockwise).
- Check the cleanliness of the shaft and the rotor.
- Place the rotor on the shaft.
- Screw the shaft nut and secure it by using the key provided in the tool kit.
- In the case of a swing-out rotor, grease the trunnions.
- Place the buckets on the trunnions after checking the weight of the buckets is identical (the weight is engraved on the bottom of each bucket).

2.9 - BALANCING THE LOAD

When using fixed angle rotor, weight deviation between the tubes must not exceed ± 5 grams.


When using a swing-out rotor, the samples are placed in plastic inserts.

The load must be statically and dynamically balanced :

Static balancing consists of balancing the weight of the 2 inserts which are in opposed location in the rotor (± 25 grams).

Dynamic balancing consists of placing the tubes in each insert symmetrically according to the shaft axis of the rotor.

The following pictures give examples of correct and incorrect loading :

CAUTION : Never use old accessories on a new centrifuge.
Only use plastic ware or glassware designed for centrifugation.

CHAPTER 3

SPECIFICATIONS

3.1 - DIMENSIONS AND WEIGHT

	C3.12	CR3.12
* Height (lid open) mm/in	850/33.5	860/34
* Height (lid closed) mm/in	400/15.8	410/16.2
* Width mm/in	480/18.9	765/30.1
* Depth mm/in	580/22.8	580/22.8
* Weight kg/lb	72/158	92/202

3.2 - ELECTRICAL SPECIFICATIONS

* Mains voltage	230 V $\pm 10\%$ - 50Hz 120 V $\pm 10\%$ - 60Hz	ref. 11175085 ref. 11175087	ref. 11175088 ref. 11175090
* Motor	220 V 120 V	25509570 25509571	25509570 25509571
* Power consumption	220 V 120 V	3 A 6 A	4 A 8 A

3.3 - PERFORMANCE AND ACCESSORIES

SPEED CONTROL

Range :	0 to 10 000 r.p.m
Increments :	250 r.p.m
Braking system :	inertial or electrical
Max. braking :	1 min 30 sec from 3600 r.p.m to 0 (T4 swing-out rotor)

SPEED READOUT

Display :	0 to 9 999 r.p.m
Increments :	10 r.p.m
Accuracy :	± 100 r.p.m

ENERGY

Max. density :	1.2 g/ml
Max. load on drive shaft :	7.7 kg
Max. energy :	20380 J

The maximum rotational energy is not systematically that obtained with the maximum permissible load on the drive shaft.

TIMER

Range :	1 to 30 min + hold position
Increments :	1 min
Accuracy :	$\pm 2 \%$

TEMPERATURE CONTROL

Range :	0 to 40°C
Increments :	1°C

TEMPERATURE READOUT

Range :	0 to 50°C
Increments :	1°C
Accuracy :	$\pm 2^\circ\text{C}$

The above performance specifications are reliable for an input voltage of :

220 Vac $\pm 10\%$ - 50 Hz or
120 Vac $\pm 10\%$ - 60 Hz

ACCESSORIES

C3.12 and CR3.12 centrifuges can be equipped with the following accessories :

- 1 swing-out rotor (T4 Cat. Nbr. : 11175366) which allows the centrifugation of many types of tubes in plastic inserts fitting standard or sealed buckets.
- 6 fixed angle rotors for centrifugation at higher speed.

SWING-OUT ROTOR

FIXED-ANGLE ROTOR

The following tables give the specifications of the accessories and their performance. Another table gives the available JOUAN glassware and the insert in which they can be used.

T4 SWING-OUT ROTOR : 11175366

Bucket or carrier (set of 4)	Rack or insert (set of 4)	Rack or insert colour	Rack or insert capacity	Max. tube Ø (mm)	Min/max. tube height (mm)	Tube volume (ml)	Max. radius (mm)	C3.12		CR3.12	
								Max. rpm	Max. RCF(g)	Max. rpm	Max. RCF(g)
Standard buckets 11175367	11175340	Black	42	8.5	35/-	0.5-0.8	174	3500	2380	3600	2520
	11175341	Salmon	20	11	35/-	1.5-2.2	174	3500	2380	3600	2520
	11175342	Beige	30	11	65/107	3/5	174	3500	2380	3600	2520
	11175343	Grey	26	12.7	67/107	5	174	3500	2380	3600	2520
	11175357	White	7	18	35/-	5 flat bottom	174	3500	2380	3600	2520
	11175344	Light blue	20	13.5	65/107	6	174	3500	2380	3600	2520
	11175345	Yellow	18	13.5	65/108	7 vac.	174	3500	2380	3600	2520
Transparent sealed lid 11175368	11175346	Red	12	16.5	65/112	10 vac.	174	3500	2380	3600	2520
	11175347	Cream	12	17.5	65/112	10 flat bottom	174	3500	2380	3600	2520
	11175358	White	7	18	90/122	12 flat bottom	174	3500	2380	3600	2520
	11175348	White	7	18	90/122	15 conical	176	3500	2410	3600	2550
	11175349	Orange	7	16.5	90/127	15 vac.	174	3500	2380	3600	2520
	11175350	Brown	4	16.5	90/133	15 vac.	176	3500	2410	3600	2550
	11175351	Coral	8	22.2	90/116	20	176	3500	2410	3600	2550
Spare seal 25439074	11175352	Light yellow	5	25	90/120	25	176	3500	2410	3600	2550
	11175375	Light yellow	5	25	90/120	25 flat bottom	176	3500	2410	3600	2550
	11175374	Dark red	4	30	90/122	25/50 flat bottom	176	3500	2410	3600	2550
	11175353	Dark red	4	30.5	90/122	50	176	3500	2410	3600	2550
	11175354	Light blue	3	30	90/117	50 conical	176	3500	2410	3600	2550
	11175376	Light blue	3	30	90/117	50 flat bottom	176	3500	2410	3600	2550
	11175361	Dark blue	2	35	90/122	50	176	3500	2410	3600	2550
	11175355	Green	1	45	90/132	100	176	3500	2410	3600	2550
	11175356	Turquoise	1	63	65/125	250 flat bottom	172	3500	2350	3700	2490
11175359 (set of 2)	Black	1 cyto-bucket (11172499 : set of 2 buckets)					181	1570	500	1570	500
11175360	5 sets of 4 decanting rubber pads (Ø 10.5 - 11 - 11.5 - 11.9 - 12.1 mm) for 11175343										
Microtitration 11175377			3 plates			96 wells	158	3400	2045	3400	2045

ANGLE ROTORS

Rotor	Catalogue number	Capacity	Holes		Min./max. tube ht (mm)	Max. radius (mm)	Max. r.p.m.	Max. RCF (g)
			Diameter (mm)	Depth (mm)				
CD1	11175181	8 x 100 ml	44.5	101	106/115	134	5900	5215
CD2	11175183	10 x 50 ml	35	92	97/115	130	6000	5232
CD3	11175186	40 x 15 ml	18.5	93	96/110 96/110	132 145	5700	4795 5267
CD4	11175180	40 x 15 ml	16.5	93	96/131 97/131	131 144	5700	4758 5231
CD5 Sealed	11175187	40 x 15 ml	16.5	110 107	115/130 112/127	136 151	5600	4768 5294
Eppendorf	11174217	30 x 1.5/2 ml	11	open	40	115	6500	5432

JOUAN GLASSWARE

Capacity	Material	Bottom	Type	Standard tube cat. number	Graduated tube cat. number	Insert cat. number	Rotor model
12 ml	GLASS	conical	O	41191314	41191315	11175351	
12 ml	GLASS	round	C	41191212		11175351	
16 ml	GLASS	round	O	41191110	41192111	11175351	
24 ml	PP	round	C	41193113		11175351	
20 ml	GLASS	conical	C	41191238		11175361	
27 ml	GLASS	conical	O	41191336	41192337	11175353	CD2
40 ml	GLASS	round	C	41191232		11175361	CD2
47 ml	GLASS	round	O	41191130	41192131	11175361	CD2
58 ml	PP	round	C	41193134		11175361	CD2
66 ml	GLASS	conical	C	41191466		11175355	CD1
95 ml	GLASS	round	C	41191262		11175355	CD1
104 ml	GLASS	round	O	41191160		11175355	CD1
130 ml	PP	round	C	41193164		11175355	CD1

PP : Polypropylene O : Open C : Capped

CHAPTER 4

OPERATING PRINCIPLES

4.1 - DRIVE SYSTEM

The mains input voltage is converted to a variable DC voltage. This DC voltage is varied by comparison of the speed control voltage setting with the tachometer feedback voltage. The greater the DC voltage applied to the motor, the faster the motor will spin a given rotor.

The C3.12 and CR3.12 centrifuges feature current limiting dependent on motor load by rotor and accessories.

The shape of the rotor is primary load factor : a horizontal rotor presents a greater load than a fixed-angle rotor. The larger the diameter of a rotor, with accessories, the greater the load, and the lower the maximum R.P.M.

Relative centrifugal force (RCF), at the circumference of a rotor and bucket combination, is directly proportional to the (R.P.M.) and radius of the rotor.

Therefore, a greater R.P.M. and/or a larger radius produces a greater R.C.F. and improved separation of substances.

4.2 - SAFETY INTERLOCK SYSTEM

The C3.12 and CR3.12 are equipped with an interlock system that prevents opening of the centrifuge lid when the rotor is spinning. The centrifuge will not operate until the lid is closed and latched in place.

The lid remains latched until the rotor stops spinning. The lid access lamp indicates when the rotor has stopped and then the button on the right hand side of the centrifuge and the LID key can be actuated to open the lid.

NOTE : If a power failure occurs, access to the samples in the centrifuge is possible. Insert and push the rod provided with the centrifuge accessory kit, into the access hole located on the top left side of the centrifuge. Perform bypass only under emergency conditions.

4.3 - IMBALANCE DETECTION SYSTEM

A red light located on the front panel indicates when an excessive imbalance has been detected. The centrifuge is then in braking phase. Open the lid for check the load. To deactivate the imbalance detection light, switch OFF the centrifuge.

As long as the light is ON, it is impossible to start the motor.

4.4 - REFRIGERATION SYSTEM

The CR3.12 is the refrigerated version of the C3.12.


Refrigeration is obtained by evaporating a refrigerant (0% CFC) in a copper coil wrapped around the rotor chamber. The sealed compressor is all ON/all OFF controlled.

The cooling group condenser is cooled by air ventilation provided by a fan.

CHAPTER 5

INSTRUCTIONS FOR USE

5.1 - CONTROLS AND INDICATORS

All controls and indicators are located on the front panel, except the mains switch which is on the left hand side of the rear panel, and the lid interlock button which is situated on the right hand side of the machine.

CONTROL INDICATOR	FUNCTION
"START" key (1)	Starts the centrifugation (if the timer is not set at 0) when the lid is properly locked.
"STOP" key (2) "BRAKE" key (3) and light (4)	Starts the braking phase. It is a manual stop key which is necessary to activate when using the timer on the "∞" position. Two braking rates are available : Soft braking when the light is OFF and reinforced braking when the light is ON. The key is used to select the required rate.
"LID" key (5), and light (6)	The light (6) indicates the rotor has stopped. Then the lid can be opened by pressing the LID key (5) while holding the lid interlock button (13) pressed in.
Imbalance detection light (7)	When excessive vibrations occur, imbalance is detected : automatic shutdown is activated and the red imbalance detection light comes ON. Opening and closing the lid will switch off the light.
Temperature control light (8)	On the CR3.12 only, this light indicates the cooling group compressor is energized.
Speed setting knob (9)	Speed control potentiometer : from 0 to 10000 rpm.
Digital display (10)	It indicates the actual speed with an accuracy of ± 100 at 6000 RPM.

Timer (11)	30 mn electromechanical timer. An "infinite" position allows longer runs.
Temperature setting knob (12)	It allows the rotor chamber temperature control within the range 0-40°C. The minimum temperature obtained at 3600 RPM with the T4 rotor for a continuous run is 4°C when the room temperature is 25°C.
Lid interlock button (13)	Push and hold, then press "LID" (5) key to open lid.

5.2 - PRELIMINARY OPERATIONS

- Switch on the unit.
- Open the lid.
- Install the accessories.
- Close the lid.
- Select the speed.
- Select the temperature (on CR3.12 only).
- Select the braking rate.
- Select the duration.

5.3 - CENTRIFUGATION CYCLE

The centrifuge is ready to start. Press the START key. The motor starts and the timer counts down.

The rotor reaches the set speed (1 mn from 0 to 3300 RPM with a loaded T4 rotor).

If the displayed speed does not corresponds to the desired setting, adjust it by using the speed setting knob. **Never choose a speed higher than the maximum permitted speed of each rotor** (see table page 3-3).

On the CR3.12, the cooling group running is indicated by a light.

5.4 - END OF THE CYCLE

When the timer reaches 0 or when the STOP key is activated, the centrifuge starts the braking phase, according to the selected braking rate which can be modified at any time.

When the speed becomes lower than 6 r.p.m, the "lid" light is ON and the lid can be opened. Wait until speed display shows "0 rpm" before pressing and holding the lid interlock button and then pressing the lid key.

In case of a power cut-out, the centrifuge stops by coasting. It is possible to manually open the lid in order to get back the samples (see § 2.7).

CAUTION : As the rotor may be rotating wait for its stop before opening.

CHAPTER 6

HAZARDS, PRECAUTIONS AND LIMITATIONS OF USE

6.1 - SPEED CONTROL

A view port located on the lid allow the measurement of the actual speed by a phototachometer.

Every 4 to 6 months, check for correct rpm readout and speed control setting which should be within 100 rpm of the actual.

6.2 - CAUTIONS

- Never try to by pass the lid lock safety while the rotor is spinning,
- Do not try to open the lid until speed display returns to "0 rpm",
- Only use correctly grounded mains source.

Special attention is required to the following :

- Installation of the unit : proper ventilation, horizontality of the centrifuge,
- Accessory handling : rotor buckets,
- Rotor installation : check the rotor is locked in position before use,
- Cleaning of the accessories, of the rotor chamber, lubrication of the trunnions.
- Centrifuge ware : shape, material and quality of the vessels must be in accordance with the performance of the centrifuge,
- Samples : cleaning of the accessories is particularly necessary when using corrosive products in the samples (saline, acids, bases),
- Load balancing.

6.3 - IEC 1010-2-020 STANDARD

The IEC 1010-2-020 standard relative to laboratory centrifuges recommends the user :

- To mark out a clearance envelope 30 cm around the centrifuge or establish special procedures for the non entry of all persons or all dangerous materials into this space during the operation of the centrifuge.
- To provide an emergency stop switch enabling the mains power to be cut in the case of a malfunction. This switch should be placed at a safe distance from the centrifuge, preferably in a different room from that in which the centrifuge is situated.

6.4 - OPERATIONAL LIMITATIONS

CAUTION : C3.12 - CR3.12 CENTRIFUGES ARE NOT EXPLOSION PROOF.

6.5 - SAMPLES

In order to get a good separation of the substances, the user must refer to chapter 1 to set the parameter values according to the vessels and to the accessory used.

A soft braking may be required to avoid resuspension of the sample.

6.6 - HAZARDS

Never service the centrifuge without disconnecting the mains.

High voltage is present behind the front panel.

6.7 - IMPROPER USE


- Use only rotors and accessories designed for use in the C3.12 and CR3.12.
- Do not attempt to override the lid interlock assembly.
- Load the rotor only in the recommended way (chapter 2). As the centrifuge starts to spin, an improperly loaded rotor could cause enough force to damage the rotor chamber.

6.8 - BIOLOGICAL RISKS

By the action of the turbulence created in the bowl a centrifuge is an aerosol generator. In the case of products presenting a biological risk you must take suitable precautions to prevent or reduce this risk and also train operators in the use of this instrument and its accessories in the context of their applications.

The following information should be given to all users including those persons using the equipment infrequently :

- A ventilated centrifuge (without built-in refrigeration) permanently exhausts into the room air which has passed through its bowl while a refrigerated (or thermostated) centrifuge contains a relatively sealed bowl. Thus the risks from the spread of aerosols are minimised in a refrigerated centrifuge and significant in a ventilated centrifuge.
- Aerosols are produced naturally when using uncapped tubes or bottles. They are also produced when a container deforms or breaks provoking a smaller or greater sample leak. We draw your attention to the use of tubes used in angle rotors. Even if the liquid level does not touch the cap when in a static state, it could do so during centrifugation and leak. You should therefore follow the instructions of the tube supplier concerning the fill level and sealing of the tube.

Once a biological risk is known or suspected, samples should be placed in a sealed container. Should such containers not be available, sealed accessories should be used such as buckets or angle rotors with sealed lids.

- Duration of aerosol presence : when aerosols are created in a centrifuge bowl it persists in the air for a period of 3-5 minutes after the rotor stops moving. Indeed, the action of opening the lid itself provokes the dispersion of aerosols in the environment.
- When an elevated risk of aerosols or of breakage is perceived the accessories should be handled using gloves and opened in an environment ensuring the protection of the operator (safety cabinet, glove chamber, wearing a mask ...), even if they are apparently sealed.

CHAPTER 7

SERVICING AND PREVENTATIVE MAINTENANCE

CAUTION : All cleaning should be performed with the centrifuge disconnected from power supply.

7.1 - CLEANING

Daily :

- No daily cleaning is required, except for accidental glass breakage, or a large amount of spillage in the rotor chamber.

Weekly :

- Clean the rotor chamber and the accessories with a cloth or cotton wool pad dipped in alcohol (70°). The centrifuges are equipped with a drain hose, attached to the chamber to allow flushing of the chamber if necessary.
- After cleaning the accessories, always rinse them with clear water and dry.

7.2 - MAINTENANCE

Speed readout :


Speed readout control using a photo-tachometer through the lid viewport.

Speed control :

Check that the centrifuge can attain the maximum r.p.m. possible for your rotor and accessories, as specified in chapter 1.

Motor brush replacement :


- Disconnect the unit from the mains.
- Remove the front panel.
- Unscrew the plastic brush caps and remove the brushes.
- Check the length of carbon brushes. If shorter than 13 mm (1/2 in) then replace the two brushes. If the motor brushes are in good condition, they should be relocated in their original position in the brush holders.
- Run for 2 hours at 2000 rpm after renewing the brushes in order to seat them.

7.3 - ROTOR REMOVAL

To remove the rotor from the motor shaft :

- Using the wrench supplied, remove the rotor lock nut by turning the wrench clockwise.
- Using the puller supplied, position its screw above the motor shaft end and position the puller arms to catch beneath two rotor arms.
- With the wrench supplied, tighten the puller screw clockwise until the rotor loosens on the motor shaft.

7.4 - FUSE REPLACEMENT

Changing of fuses must not be carried out by the operator but by a qualified engineer who will diagnose the fault before replacing the fuses.

7.5 - SPARE PART LIST

QUANTITY	DESCRIPTION	CAT. NUMBER
2	Carbon brush 6 x 16 mm	89000214 <input type="checkbox"/>
2	Carbon brush 6 x 10 mm	89000056 <input type="checkbox"/>
1	Fuse T 6.3 x 32 - 1 A (120 V)	86001383
1	Fuse T 6.3 x 32 - 0.5 A (230 V)	86001382

NOTE :

Any repair performed by the user, with the exception of brush replacement, cancels the warranty.