PowerWaveTM

Operator's Manual





PowerWaveTM

Microplate Spectrophotometer Operator's Manual

June 2011 © 2011 Part Number 7281000 Revision F BioTek® Instruments, Inc.

Notices

BioTek® Instruments, Inc.

Highland Park, P.O. Box 998 Winooski, Vermont 05404-0998 USA

All Rights Reserved

© 2011, BioTek® Instruments, Incorporated. No part of this publication may be reproduced, transcribed, or transmitted in any form, or by any means electronic or mechanical, including photocopying and recording, for any purpose other than the purchaser's use without written permission of BioTek Instruments, Inc.

Trademarks

BioTek® is a registered trademark, and PowerWaveTM, Gen5TM, 4-ZoneTM and BioStackTM are trademarks of BioTek Instruments, Inc. BioCellTM is a trademark of BioTek Instruments and is patented under U.S. patent number 5,963,318.

Microsoft®, and Windows®, are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries.

All other trademarks are the property of their respective holders.

Restrictions and Liabilities

Information in this document is subject to change and does not represent a commitment by BioTek Instruments, Inc. Changes made to the information in this document will be incorporated in new editions of the publication. No responsibility is assumed by BioTek for the use or reliability of software or equipment that is not supplied by BioTek or its affiliated dealers.

Contents

Preface	2
Notices	2
All Rights Reserved	2
Trademarks	2
Restrictions and Liabilities	2
Contents	3
Contact Information	5
Customer Service and Sales	5
Service/TAC	5
European Coordination Center/Authorized European Representation	esentative5
Revision History	6
Document Conventions	8
Intended Use Statement	8
Quality Control	9
Warranty and Product Registration	9
Repackaging and Shipping	9
Warnings	9
Hazards	10
Precautions	11
CE Mark	12
Directive 2004/108/EC: Electromagnetic Compatibility	12
Directive 73/23/EEC Low Voltage (Safety)	12
Directive 2002/96/EC: Waste Electrical and Electronic Equ	uipment 12
Directive 98/79/EC: In Vitro Diagnostics	13
Electromagnetic Interference and Susceptibility	13
USA FCC CLASS A	13
Canadian Department of Communications Class A	
User Safety	
Safety Symbols	15
Chaper 1: Introduction	17
Product Description	
Package Contents	18
Optional Accessories	19
Specifications	20
Microplates	20
Speed of Reading	20
Optical Specifications	21
Optical Performance	
Hardware and Environmental Specifications	22
Product Support & Service	
Contacting the Technical Assistance Center	23
Returning Instruments for Service/Repair	23

Chapter 2: Installation	25
Product Registration	26
1: Unpack and Inspect the Instrument	26
2: Remove the Carrier Shipping Bracket	26
3: Select an Appropriate Location	27
4: Connect the Power Supply	28
5: Connect the Host Computer	29
6: Install Gen5	29
7: Turn on the Reader	29
Operational/Performance Qualification	30
Repackaging and Shipping Instructions	31
Chapter 3: Operation	35
Operating the PowerWave	
Getting Started with Gen5	36
Recommendations for Optimum Performance	37
Where to Go Next	
Chapter 4: Instrument Qualification	
Overview	40
Recommended Qualification Schedule	
Qualification Procedures	
System Test	
Sample System Test Report	
Absorbance Plate Test	
Test Plate and Certificates	
Enter the Absorbance Test Plate Data	
Run the Absorbance Plate Test	
Sample Absorbance Plate Test Report	
Test Results	
Liquid Testing	
Absorbance Liquid Test 1	
Absorbance Liquid Test 2	
Absorbance Liquid Test 3 (Optional)	
Chapter 5: Decontamination	
Purpose	
Clean Plate Carrier and Exposed Surfaces	
Decontamination	
Tools and Supplies	
Chapter 6: Troubleshooting	
Error Codes	
Error Codes During Operation with the BioStack	
General Errors	
Fatal Errors	
Chapter 7: Instrument Dimensions	
Instrument Dimensions	70

Contact Information

BioTek® Instruments, Inc.

Highland Park, P.O. Box 998

Winooski, Vermont 05404-0998 USA

Customer Service and Sales

Internet: www.biotek.com

Phone: 888-451-5171 (toll free in the U.S.)

802-655-4740 (outside the U.S.)

Fax: 802-655-7941

E-Mail: customercare@biotek.com

Service/TAC

Phone: 800-242-4685 (toll free in the U.S.)

802-655-4740 (outside the U.S.)

802-654-0638 Fax:

E-Mail: tac@biotek.com

European Coordination Center/Authorized European Representative

BioTek® Instruments GmbH

Kocherwaldstrasse 34

D-74177 Bad Friedrichshall

Germany

Internet: www.biotek.de

Phone: +49 (0) 7136 9680

Fax: +49 (0) 7136 968 111

E-Mail: info@biotek.de

Revision History

Revision	Date	Changes
А	12/2001	First issue.
В	2/2002	For the Optical Performance Specifications, changed the maximum allowable Gain on Optics from 6.0 to 10.0 (page 1-5). Replaced Figures 3-1a and 3-1b, Sample Output for the System Test (pages 3-4 and 3-5). Updated technical support contact information (pages iii, 1-6, and 1-7).
		Made editorial changes.
С	8/2003	Updated TAC information in Preface and Ch. 1, and Customer Service and European contact information, Hazards, Standards, Safety Symbols, Intended Use Statement, and Warranty in Preface. Updated Introduction, Variations, Hardware Features, Software Features,
		Package Contents, Optional Accessories, and Specifications, Ch. 1. Enhanced Unpacking/Repackaging sections, Ch. 2; added Fig. 2-5, 2-6.
		Revised Ch. 3 sections: System Test, Universal Plate Test, Liquid Testing.
		Enhanced Decontamination section; added Cleaning section in App. A.
		Enhanced section on KC4 [™] , added section on KCjunior [™] in App. B, and replaced previous serial control section with a reference to Bio-Tek's serial communication protocol specification 7266201-SP.
		Replaced Error Codes tables in App. C with more detailed tables from the PowerWave™ HT Service Manual.
		Added information about the optional barcode scanner in Ch. 1, and added new App. E, Barcode Scanner.
		Referenced new versions of PowerWave™, PowerWave™ 340, and PowerWave™ HT 340 throughout manual.
		Included references to KCjunior as an additional primary operating software and substitute for KC4.
		Included information about the compatibility of the PowerWave with the Bio-Stack™ Automated Microplate Stacking System, and brief paragraphs concerning installation, serial cable connections, error codes, and alignment when operating the PowerWave with the Bio-Stack; referenced the Bio Stack Operator's Manual.
D	5/2006	Updated primarily to support introduction of Gen5™ Software.
		General: Added Gen5 references and instructions wherever KC4™ and KCjunior™ references and instructions were present. Changed 'Bio-Tek' to 'Bio-Tek,' 'Bio-Stack™ Automated Microplate Stacking System' to 'Bio-Stack Microplate Stacker,'and 'Abs' to 'OD'. Removed 'Scanning' from 'PowerWave™ Microplate Scanning Spectrophotometer'. Cover: Replaced existing cover with new design.

Revision	Date	Changes
(D)		Preface: Updated contact information, Warnings, Hazards, Pre-cautions, Safety Symbols. Removed Warranty and Registration Card. Chapter 1, Introduction: Added clarification (in Internal Barcode Scanner section) that some older models of the reader may include the scanner. Updated Package Contents, Optional Accessories, and replaced previous Technical Support pages with a Product Support and Service page. Chapter 2, Installation: Rearranged installation steps to better reflect actual practice. Chapter 3, Performance Verification/Qualification Tests: Added Gen5™ instructions for the Self Test and Absorbance Plate Test. In Recommended Qualification Schedule, moved Absorbance Plate Test and Liquid Tests from IQ to initial/annual OQ, changed PQ semiannual frequency to quarterly, and clarified criteria for running Liquid Tests 1, 2, or 3. Changed 'Universal' to 'Absorbance' in 'Universal Test Plate' and 'Universal Plate Test'. In Liquid Test 1, added BioTek wetting agent (7773002) to list of ingredients. In Liquid Test 3, changed Sigma® 'P 3563 packets' to 'PBS tablets (#4417, or equivalent).' In Liquid Tests 1 and 3, changed 'Analytical balance' to 'Precision balance.' In all Liquid Tests, added 'Corning" to "Costar' (microplates), and added note to shake plate or wait after pipetting and before reading the plate. Appendix A, Decontamination and Cleaning: Corrected dilution mixtures for bleach on page A-3 by changing '20:1' ratio for commercial bleach to '1:20', and '10:1' ratio for household bleach to '1:10'. Appendix B, Computer Control: Added new section, "Controlling the Reader with Gen5."
E	12/2009	Throughout: Removed references to models 'PowerWave' and 'PowerWave HT 340' (PowerWave HT and PowerWave 340 remain). Emphasized use of Gen5 instead of KC4 and KCjunior (which are no longer available from BioTek). Removed references to the ActiveX component. Preface: Updated Trademarks, Intended Use Statement, Hazards, Precautions, CE Mark information, and Safety Symbols. Removed lists of illustrations and tables. Ch 1 Introduction: Removed 'Variations' and 'Internal Barcode Scanner'. Updated Package Contents and Optional Accessories. Ch 2 Installation: Simplified unpacking and setup instructions. Removed Serial Pinout Description. Ch 3 Operation: New chapter. Ch 4 Instrument Qualification: Moved recommendations for optimum performance to new chapter 3. Clarified instructions for the various qualification tasks. Former Appendix B, Computer Control: Deleted this section. Moved Gen5 instructions to new chapter 3. Former Appendices A and C: Changed to Chapters 5 and 6. Former Appendices D and E: Changed to Appendices A and B.
F	5/2011	General: Removed references to outdated software KC4 and KCjunior. Updated Gen5 instructions for Gen5 version 2.x. Introduction: Deleted redundant "Hardware Features" and "Software Features" sections. Liquid Testing: Updated Liquid Test 3; removed instructions for creating the 10x concentration PBS solution. Appendices: Removed former Appendix B: Barcode Scanner.

Document Conventions

This manual uses the following typographic conventions:

\triangle	This icon calls attention to important <i>safety</i> notes.
Warning!	A <i>Warning</i> indicates the potential for bodily harm and tells you how to avoid the problem.
Caution	A <i>Caution</i> indicates potential damage to the instrument and tells you how to avoid the problem.
Note:	Bold text is primarily used for emphasis.
(i)	This icon calls attention to <i>important information</i> .

Intended Use Statement

The PowerWave is an eight-channel, automated, benchtop, general-purpose Microplate Spectrophotometer that performs optical density measurements of samples in a microplate format. The user must evaluate this instrument with PC-based software in conjunction with the specific assay. This evaluation must include the confirmation that performance characteristics for the specific assay are met.

- BioTek Gen5 software package provides the user with instrument control and data reduction capabilities.
- This product can operate with standard robotic systems, such as the BioStack Microplate Stacker.
- This product may be used for In Vitro Diagnostic, research and development, or other non-clinical purposes.

Quality Control

It is considered good laboratory practice to run laboratory samples according to instructions and specific recommendations included in the package insert or standard laboratory protocol for the test to be conducted. Failure to conduct Quality Control checks could result in erroneous test data.

Warranty and Product Registration

Review the Warranty information that shipped with your product. Register your product(s) with BioTek to ensure that you receive important information and updates. Contact the Customer Resource Center (CRC) at www.biotek.com or by calling 888/451-5171 or 802/655-4740.

Repackaging and Shipping



If you need to ship the instrument to BioTek for service or repair, contact BioTek for a **Return Materials Authorization (RMA)** number, and be sure to use the original packing materials. Other forms of commercially available packaging are not recommended and can void the warranty. If the original packing materials have been damaged or lost, contact BioTek for replacement packing.

Warnings



Operate the instrument on a level surface away from excessive humidity.

Strong light can reduce the linear performance range of the instrument.

Measurement values may be affected by extraneous particles (such as dust) in the microplate wells. A clean work area is necessary to ensure accurate readings.

When operated in a safe environment according to the instructions in this document, there are no known hazards associated with the instrument. However, the operator should be aware of certain situations that could result in serious injury; these may vary depending on the instrument model. See **Hazards** and **Precautions**.

The following hazard warnings are provided to help avoid injury:



Warning! Power Rating. The instrument's power supply cord must be connected to a power receptacle that provides voltage and current within the specified rating for the system. Use of an incompatible power receptacle may produce electrical shock and fire hazards.

Warning! Electrical Grounding. Never use a two-prong plug adapter to connect primary power to the external power supply. Use of an adapter disconnects the utility ground, creating a severe shock hazard. Always connect the power supply directly to an appropriate receptacle with a functional ground.



Warning! Internal Voltage. Always turn off the power switch and unplug the power supply before cleaning the outer surface of the instrument.

Warning! Liquids. Avoid spilling liquids on the reader; fluid seepage into internal components creates a potential for shock hazard or instrument damage. If a spill occurs while a program is running, abort the program and turn off the instrument. Wipe up all spills immediately. Do not operate the instrument if internal components have been exposed to fluid.

Warning! Unspecified Use. Failure to operate this equipment according to the guidelines and safeguards specified in this manual could result in a hazardous condition.

Warning! Software Quality Control. The operator must follow the manufacturer's assay package insert when modifying software parameters and establishing reading methods. Failure to conduct quality control checks could result in erroneous test data.

Warning! Reader Data Reduction Protocol. No limits are applied to the raw absorbance data. All information exported via computer control must be thoroughly analyzed by the operator.



Warning! Potential Biohazards. Some assays or specimens may pose a biohazard. Adequate safety precautions should be taken as outlined in the assay's package insert. Always wear safety glasses and appropriate protective equipment, such as chemically resistant rubber gloves and an apron.

Precautions

The following precautions are provided to help avoid damage to the instrument:



Caution: Service. The instrument should be serviced by BioTek authorized service personnel. Only qualified technical personnel should perform troubleshooting and service procedures on internal components.

Caution: Environmental Conditions. Do not expose the system to temperature extremes. For proper operation, ambient temperatures should remain within the range listed in the Specifications section of Chapter 1. Performance may be adversely affected if temperatures fluctuate above or below this range. Storage temperature limits are broader.

Caution: Sodium Hypochlorite. Do not expose any part of the instrument to the recommended diluted sodium hypochlorite solution (bleach) for more than 20 minutes. Prolonged contact may damage the instrument surfaces. Be certain to rinse and thoroughly wipe all surfaces.

Caution: External Power Supply. Only use the power supply shipped with the instrument. Operate this power supply within the range of line voltages listed on it.

Caution: Carrier Shipping Bracket. The microplate carrier shipping bracket must be removed before operating the instrument and reinstalled before repackaging the instrument for shipment.

Caution: Disposal. This instrument contains printed circuit boards and wiring with lead solder. Dispose of the instrument according to Directive 2002/96/EC, "on waste electrical and electronic equipment (WEEE)" or local ordinances.

Caution: Warranty. Failure to follow preventive maintenance protocols may void the warranty. See **Chapter 5** for preventive maintenance procedures.

Caution: Electromagnetic Environment. Per IEC 61326-2-6 it is the user's responsibility to ensure that a compatible electromagnetic environment for this instrument is provided and maintained in order that the device will perform as intended.

Caution: Electromagnetic Compatibility. Do not use this device in close proximity to sources of strong electromagnetic radiation (e.g., unshielded intentional RF sources), because these may interfere with the proper operation.

CE Mark



Based on the programs described below and information contained herein, this product bears the CE mark.

❖ See the Declaration of Conformity for more information.

Directive 2004/108/EC: Electromagnetic Compatibility

Emissions—CLASS A

The system has been type-tested by an independent, accredited testing laboratory and found to meet the requirements of EN 61326-1 and EN 61326 2 6: Class A for Radiated Emissions and Line Conducted Emissions.

Verification of compliance was conducted to the limits and methods of EN 55011 – (CISPR 11) Class A. In a domestic environment it may cause radio interference, in which case you may need to mitigate the interference.

Immunity

The system has been type-tested by an independent, accredited testing laboratory and found to meet the requirements of EN 61326-1 and EN 61326-2-6 for Immunity. Verification of compliance was conducted to the limits and methods of the following:

EN 61000-4-2, Electrostatic Discharge

EN 61000-4-3, Radiated EM Fields

EN 61000-4-4, Electrical Fast Transient/Burst

EN 61000-4-5, Surge Immunity

EN 61000-4-6, Conducted Disturbances from RFI

EN 61000-4-11, Voltage Dips, Short Interruptions and Variations

Directive 73/23/EEC Low Voltage (Safety)

The system has been type-tested by an independent testing laboratory and was found to meet the requirements of EC Directive 73/23/EEC for Low Voltage. Verification of compliance was conducted to the limits and methods of EN 61010-1, "Safety requirement for electrical equipment for measurement, control and laboratory use. Part 1, General requirements."

Directive 2002/96/EC: Waste Electrical and Electronic Equipment

Disposal Notice: This instrument contains printed circuit boards and wiring with lead solder. Dispose of the instrument according to Directive 2002/96/EC, "on waste electrical and electronic equipment (WEEE)" or local ordinances.

Directive 98/79/EC: In Vitro Diagnostics

- Product registration with competent authorities.
- Traceability to the U.S. National Institute of Standards and Technology (NIST): Optical density measurements are traceable to NIST.

Electromagnetic Interference and Susceptibility

USA FCC CLASS A

Warning: Changes or modifications to this unit not expressly approved by the manufacturer could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules.

These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. Like all similar equipment, this equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause interference, in which case the user will be required to correct the interference at his own expense.

Canadian Department of Communications Class A

This digital apparatus does not exceed Class A limits for radio emissions from digital apparatus set out in the Radio Interference Regulations of the Canadian Department of Communications.

Le present appareil numerique n'emet pas de bruits radioelectriques depassant les limites applicables aux appareils numerique de la Class A prescrites dans le Reglement sur le brouillage radioelectrique edicte par le ministere des Communications du Canada.

User Safety

This device has been type-tested by an independent laboratory and found to meet the requirements of the following:

North America

- Canadian Standards Association CAN/CSA C22.2 No. 61010-1, "Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use, Part 1: General Requirements"
- UL 61010-1, "Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use, Part 1: General Requirements"

International

• EN 61010-1, "Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use, Part 1: General Requirements"

Safety Symbols

Some of these symbols appear on the instrument or accessories:

Alternating current Courant alternatif Wechselstrom Corriente alterna Corrente alternata	Both direct and alternating current Courant continu et courant alternatif Gleich - und Wechselstrom Corriente continua y corriente alterna Corrente continua e corrente alternata
Direct current Courant continu Gleichstrom Corriente continua Corrente continua	Earth ground terminal Borne de terre Erde (Betriebserde) Borne de tierra Terra (di funzionamento)
On (Supply) Marche (alimentation) Ein (Verbindung mit dem Netz) Conectado Chiuso	Protective conductor terminal Borne de terre de protection Schutzleiteranschluss Borne de tierra de protección Terra di protezione
Off (Supply) Arrêt (alimentation) Aus (Trennung vom Netz) Desconectado Aperto (sconnessione dalla rete di alimentazione)	Caution (refer to accompanying documents) Attention (voir documents d'accompanement) Achtung siehe Begleitpapiere Atención (vease los documentos incluidos) Attenzione, consultare la doc annessa
Warning, risk of electric shock Attention, risque de choc électrique Gefährliche elektrische schlag Precaución, riesgo de sacudida eléctrica Attenzione, rischio di scossa elettrica	Warning, risk of crushing or pinching Attention, risque d'écrasement et pincement Warnen, Gefahr des Zerquetschens und Klemmen Precaución, riesgo del machacamiento y sejeción Attenzione, rischio di schiacciare ed intrappolarsi
Warning, hot surface Attention, surface chaude Warnen, heiße Oberfläche Precaución, superficie caliente Attenzione, superficie calda	Warning, potential biohazards Attention, risques biologiques potentiels Warnung! Moegliche biologische Giftstoffe Atención, riesgos biológicos Attenzione, rischio biologico



In vitro diagnostic medical device

Dispositif médical de diagnostic in vitro Medizinisches In-Vitro-Diagnostikum Dispositivo médico de diagnóstico in vitro Dispositivo medico

diagnostico in vitro



Separate collection for electrical and electronic equipment

Les équipements électriques et électroniques font l'objet d'une collecte sélective

Getrennte Sammlung von Elektro- und Elektronikgeräten

Recogida selectiva de aparatos eléctricos y electrónicos

Raccolta separata delle apparecchiature elettriche ed elettroniche



Consult instructions for use Consulter la notice d'emploi Gebrauchsanweisung beachten Consultar las instrucciones de uso Consultare le istruzioni per

Introduction

This chapter introduces the PowerWave Microplate Spectrophotometer, describes its features and specifications, and provides contact information for technical assistance.

Product Description	18
Package Contents	18
Optional Accessories	19
Specifications	20
Product Support and Service	23

Product Description

The PowerWave offers tunable wavelength selection and wavelength scanning without the need for interference filters. The eight-channel reader is computer controlled via BioTek's Gen5 software. Key features include:

- A variety of read methods including endpoint, kinetic, multiwavelength, and spectral scanning.
- A monochromator for continuous wavelength selection from 200 to 999 nm (or 340 to 999 nm for the PowerWave 340), in 1-nm increments
- A xenon flash lamp for both UV and visible light absorbance measurements.
- Superior optical specifications, with an extended dynamic range of up to 4.000 OD.
- Ability to read standard 96- and 384-well (PowerWave HT) microplates and BioTek's patented BioCell quartz vessel for 1 cm measurements.
- Three reading speeds: normal, rapid, and sweep mode.
- A unique 4-ZoneTM temperature control from 4° over ambient to 50°C that ensures superior temperature uniformity necessary for kinetic assays.
- Low, medium, high, and variable plate shake speeds with adjustable durations.
- Robot accessible carrier. Compatible with BioTek's BioStack Microplate Stacker.
 - If you purchased the BioStack to operate with the PowerWave, refer to the BioStack Operator's Manual for installation, setup, and operation instructions. If you are interested in purchasing the BioStack, contact your local BioTek dealer or visit www.biotek.com.

Package Contents

- ❖ Package contents and part numbers are subject to change. Please contact BioTek Customer Care if you have any questions.
- Gen5 Software (PN 5320200)
- Power supply (PN 76053) and power cord (PN varies according to country of use)
- Serial cable (PN 75053)
- PowerWave Operator's Manual (PN 7281000)

Optional Accessories

- Accessories and part numbers are subject to change. Please contact BioTek Customer Care if you have any questions, or visit www.biotek.com and use the Accessories search feature.
- USB to Serial Adapter (PN 75104)
- Absorbance Test Plate (PN 7260522)
- BioCell Quartz vessel for 1 cm wavelength fixed pathlength absorbance measurements (PN 7272051); adapter plate for up to eight BioCells (PN 7270512)
- PowerWave Product Qualification (IQ-OQ-PQ) package (7280520)
- Liquid Test Solutions:
 - ➤ BioTek Wetting Agent (PN 7773002)
 - ➤ BioTek QC Check Solution No. 1 (PN 7120779 for 25 ml; PN 7120782 for 125 ml)
- BioStack Microplate Stacker (contact Customer Care)

Specifications

Microplates

- All models accommodate standard 96-well microplates, and up to 8 BioCells.
- The PowerWave HT also accommodates standard 384-well microplates.

Speed of Reading

The plate read time and accuracy are dependent on the method of reading:

- Normal mode is the slowest of the three available modes. After positioning the well over the beam, the instrument waits 100 milliseconds before taking the measurement (8-flash data collection). Note: The 100 ms delay is to allow for more complete fluid settling.
- Rapid mode is faster than Normal mode because the instrument does not wait before taking the measurement (8-flash data collection).
- Sweep is the fastest of the three available modes. The plate carrier sweeps each row past the optics channel without stopping, and collects data with a single flash at each well as it goes by.

The following read times are based on a single or dual wavelength measurement. Actual reading speeds may vary, depending upon the reading wavelength selected. Each wavelength has a unique location within the monochromator, and the different locations require varying amounts of time to position.

96-Well Read Timing	630 nm	630/450 nm
Normal Read Mode	Single	Dual
Endpoint	16 to 25 sec.	26 to 44 sec.
Rapid Read Mode	Single	Dual
Endpoint	16 sec.	26 sec.
Sweep Read Mode	Single	Dual
Endpoint	11 sec.	16 sec.

Kinetics: All three read modes are available in Kinetics mode. Single wavelength reads are limited to the following minimum times.

20 seconds from A1 to A1 in Normal mode, single wavelength, depending upon density of solution.

11 seconds from A1 to A1 in Rapid mode, single wavelength.

5 seconds from A1 to A1 in Sweep mode, single wavelength.

384-Well Read Timing	630 nm	630/450 nm
Normal Read Mode	Single	Dual
Endpoint	32 to 67 sec.	57 to 129 sec
Rapid Read Mode	Single	Dual
Endpoint	28 sec.	49 to 51 sec.
Sweep Read Mode	Single	Dual
Endpoint	17 sec.	28 sec.

Kinetics: All three read modes are available in Kinetics mode. Single wavelength reads are limited to the following minimum times.

66 seconds from A1 to A1 in Normal mode, single wavelength, depending upon density of solution.

23 seconds from A1 to A1 in Rapid mode, single wavelength.

11 seconds from A1 to A1 in Sweep mode, single wavelength.

Optical Specifications

λ range: 200 to 999 nm (PowerWave HT)

340 to 999 nm (PowerWave 340)

 λ accuracy: ± 2 nm

 λ repeatability: ± 0.2 nm

λ bandpass: 5 nm

Optical Performance

Flat- and round-bottom full-well plates

Absorbance Measurement Range: 0.000 to 4.000 OD

Accuracy:

0.000 to 2.000 OD $\pm 1.0\% \pm 0.010$ OD Normal and Rapid modes, all plates 2.000 to 2.500 OD $\pm 3.0\% \pm 0.010$ OD Normal and Rapid modes, all plates $2.500 \text{ to } 3.000 \text{ OD} \pm 3.0\% \pm 0.010 \text{ OD Normal 96-well plates only}$ 0.000 to 1.000 OD $\pm 1.0\% \pm 0.010$ OD Sweep mode, all plates

Linearity:

0.000 to 2.000 OD ± 1.0% Normal and Rapid modes, 96-well plates 0.000 to 2.000 OD $\pm 2.0\%$ Normal and Rapid modes, 384-well plates 2.000 to 2.500 OD $\pm 3.0\%$ Normal and Rapid modes, all plates 2.500 to 3.000 OD $\pm 3.0\%$ Normal mode, 96-well plates only 0.000 to 1.000 OD $\pm 1.0\%$ Sweep mode, all plates

Repeatability:

0.000 to 2.000 OD $\pm 1.0\% \pm 0.005$ OD Normal and Rapid modes, 96-well plates 0.000 to 2.000 OD $\pm 2.0\% \pm 0.010$ OD Normal and Rapid modes, 384-well plates 2.000 to 2.500 OD $\pm 3.0\% \pm 0.005$ OD Normal and Rapid modes, all plates $2.500 \text{ to } 3.000 \text{ OD} \pm 3.0\% \pm 0.005 \text{ OD Normal mode, } 96\text{-well plates only}$ $0.000 \text{ to } 1.000 \text{ OD} \pm 2.0\% \pm 0.010 \text{ OD Sweep mode, all plates}$

For the above performance, the Gain on Optics test should be below 10.0.

Hardware and Environmental Specifications

Light Source: Xenon flash light source

- 10 W max. average power

- Life: 1 billion flashes

Dimensions: 16.0" deep x 8.5" wide x 8.5" high

(40.6 cm deep x 21.6 cm wide x 21.6 cm high)

Weight: 24 lb. (10.9 kg)

Operational temperature 18° - 40°C **Environment:**

Humidity: 10% to 80%, non-condensing

Power Source: 24-volt external power supply compatible with

 $100-240 \text{ V} \sim \pm 10\% \text{ } \text{@} 50-60 \text{ Hz}$

Power Consumption: 100 VA max

Temperature Control: 4°C over ambient to 50°C

Temperature Variation: ± 0.5C across the plate @ 37°C (250 μL per well with the

plate sealed)

Product Support & Service

Contacting the Technical Assistance Center

If your instrument(s) or software fails to function properly, if you have questions about how to use or maintain your products, or if you need to send an instrument to BioTek for service or repair, please contact our Technical Assistance Center (TAC).

The TAC is open from 8:30 AM to 5:30 PM (EST), Monday through Friday, excluding standard U.S. holidays.

E-Mail: tac@biotek.com **Phone:** (800) 242-4685 or **Fax**: (802) 654-0638 Web: www.biotek.com (802) 655-4740

Please be prepared to provide the following information:

- Your name and company information, along with a daytime phone or fax number, and/or an e-mail address
- The product name, model, and serial number
- The onboard software part number and version (available via Gen5 by selecting System > Reader Control > Information)
- Gen5 software version information (Help > About Gen5)
- For troubleshooting assistance or instruments needing repair, the specific steps that produce your problem and any error codes displayed in Gen5 (see also Chapter 6)

Returning Instruments for Service/Repair

If you need to return an instrument to BioTek for service or repair, please contact the TAC for a Return Materials Authorization (RMA) number and the shipping address. Repackage the instrument properly (see Chapter 2), write the RMA number on the shipping box, and ship to BioTek.

Chapter 2 Installation

This chapter provides instructions for unpacking and setting up the PowerWave.

Product Registration	26
1: Unpack and Inspect the Instrument	
·	
2: Remove the Carrier Shipping Bracket	
3: Select an Appropriate Location	27
4: Connect the Power Supply	28
5: Connect the Host Computer	29
6: Install Gen5	29
7: Turn on the Reader	29
Operational/Performance Qualification	30
Repackaging and Shipping Instructions	31

Product Registration

Register your product(s) with BioTek to ensure that you receive important information and updates. Register online through BioTek's Customer Resource Center (CRC) at www.biotek.com or by contacting BioTek Customer Care.

1: Unpack and Inspect the Instrument



Important! Save all packing materials. If you are sending the reader to BioTek for repair or replacement, use the original packing materials. Using other forms of commercially available packing materials, or failure to follow the packaging instructions at the end of this chapter, may void your warranty. If the original materials have been damaged or lost, replacements are available from BioTek (PN 7283000).

Carefully unpack the reader and accessories. Retain the packing materials for future use.

Inspect the shipping box(es), reader, and accessories for signs of damage.

If the reader is damaged, notify the carrier and your manufacturer's representative. Keep the shipping cartons and packing material for the carrier's inspection. The manufacturer will arrange for repair or replacement of your reader.

See **Preparing the PowerWave for Shipment** at the end of this chapter for repacking and shipping instructions.

2: Remove the Carrier Shipping Bracket



Important! The PowerWave ships with a microplate carrier shipping bracket that must be removed before the reader is used. See **Figure 1**.

- 1. On the front of the reader, pull down the door to the carrier compartment.
- 2. Using a screwdriver, remove the three screws that secure the shipping bracket.
- 3. Secure the bracket to the back of the reader for storage.




Figure 1: Remove the carrier shipping bracket and store it on the back of the reader

3: Select an Appropriate Location

Install the reader on a level surface in an area where ambient temperatures remain between 18°C (64°F) and 40°C (104°F). The reader is sensitive to extreme environmental conditions; avoid these conditions:

- **Excessive humidity**: Condensation directly on the sensitive electronic circuitry can cause the reader to fail internal self-checks.
- **Excessive ambient light:** Strong light can reduce the linear performance range of the reader.
- **Dust**: Optical density readings may be affected by extraneous particles in the microplate wells. A clean work area is necessary to ensure accurate readings.

4: Connect the Power Supply



Warning! Power Rating. The power supply must be connected to a power receptacle that provides voltage and current within the specified rating for the system. Use of an incompatible power receptacle may produce electrical shock and fire hazards.

Warning! Electrical Grounding. Never use a plug adapter to connect primary power to the power supply. Use of an adapter disconnects the utility ground, creating a severe shock hazard. Always connect the power supply cord directly to an appropriate receptacle with a functional ground.

- Connect the power cord to the power supply.
- Locate the power inlet on the back of the reader. Insert the power supply's plug into the reader's inlet. Tighten the plug barrel.
- Plug the cord into an appropriate power receptacle.




Figure 2: Serial port and power supply inlet

5: Connect the Host Computer

- 1. Turn the computer off. If the reader is on, turn it off.
- 2. Connect one end of the supplied serial cable to an appropriate port on the computer.
 - BioTek offers a USB to Serial Adapter (PN 75104) to connect the serial cable to a USB port on your computer. Contact BioTek Customer Care.
- 3. Connect the other end of the cable to the serial port on the back of the reader.

6: Install Gen5



The PowerWave is controlled by Gen5 software running on a host computer. There is a certain sequence of events that *must* be followed to ensure that the software is properly installed and configured. Please follow the instructions provided in *Gen5 Getting Started Guide* to install the software.

7: Turn on the Reader

- 1. Locate the power on/off switch on the front of the instrument, below the carrier eject button. See **Figure 3** on the next page. Turn on the power. The reader will perform an internal self-test and carrier homing sequence.
- 2. Verify that the following occur while the reader performs the self-test:
 - The carrier should eject outside the PowerWave, then retract to its home position inside the reader before it ejects again.
 - The LED light on the switch should remain illuminated while the power is
- 3. Press the carrier eject button. The carrier should retract and the door should close. Press it again; the carrier should eject.
- 4. If the test completes successfully, the reader is ready for use.
- 5. If the test fails, the reader will "beep" continuously. Press the carrier eject button to stop the beeping. Run a System Test using Gen5 to retrieve an error code (see Chapter 6). If the test continues to fail, contact BioTek.




Figure 3: Power on/off switch and carrier eject button

Operational/Performance Qualification

Your PowerWave Microplate Spectrophotometer was fully tested at BioTek prior to shipment and should operate properly following the successful completion of the installation and setup procedures described throughout this chapter.

If you suspect that problems occurred during shipment, if you received the reader back from BioTek following service or repair, and/or if regulatory requirements dictate that Operational/Performance Qualification is necessary, turn to *Chapter 4*, *Instrument* Qualification now to learn about BioTek's recommended OQ/PQ procedures for the PowerWave.

An Installation/Operational/Performance Qualification (IQ/OQ/PQ) package for the PowerWave is available for purchase (PN 7280520). Contact your local BioTek dealer for more information.

Repackaging and Shipping Instructions

Important! Please read all of the information provided below before preparing the PowerWave for shipment.



If the reader has been exposed to potentially hazardous material, decontaminate it to minimize the risk to all who come in contact with the reader during shipping, handling, and servicing. Decontamination prior to shipping is required by the U.S. Department of Transportation regulations. See *Chapter 5* for decontamination instructions.

Remove the microplate from the carrier before shipment. Spilled fluids can contaminate the optics and damage the instrument.



The instrument's packaging design is subject to change. If the instructions in this section do not apply to the packaging materials you are using, please contact BioTek's Technical Assistance Center for guidance.

Replace the shipping hardware before repackaging the reader. Please contact BioTek if you have misplaced the microplate carrier shipping bracket (PN 7282014) or mounting screws (3, PN 19186).

If you need to ship the reader to BioTek for service or repair, be sure to use the original packaging materials. Other forms of commercially available packaging are not recommended and can void the warranty.

The shipping materials are designed to be used no more than five times. If the original materials have been damaged, lost, or used more than five times, contact BioTek to order replacements.

- 1. Contact BioTek's Technical Assistance Center for an **RMA** (Return Materials Authorization) number and the shipping address before returning equipment for service.
- 2. Decontaminate the reader according to the instructions provided in *Chapter 4*.
- 3. Install the carrier shipping bracket (refer to *Remove the Carrier Shipping* **Bracket** on page 26):
 - If the carrier is extended, press the carrier eject button to retract it.
 - Turn off the reader.
 - Disconnect the power supply and serial cable from the back of the reader.

- Refer to **Figure 1** on page 27. Using a screwdriver, remove the carrier shipping bracket and screws from the back of the reader.
- Pull down the door to the carrier compartment.
- Install the carrier shipping bracket to the front of the carrier and mounting block.

Refer to **Figure 4** on page 33 when performing these steps:

- Place two foam caps into the bottom of the shipping container.
- 2. Slide the accessories box into the shipping container.
- Place the reader inside its plastic bag and carefully lower it into the two foam caps in the bottom of the box. Note the orientation of the reader in the box.
- 4. Place two foam caps over the reader.
- 5. Place the power supply, cord, and communication cables into the accessories box.
- Close the top of the box and secure it with shipping tape. When finished, write the RMA number on the outside of the box and ship the box to BioTek.




Figure 4: Packing the PowerWave

Chapter 3

Operation

This chapter briefly describes how to use BioTek's Gen5 software to operate the PowerWave. It also contains recommendations for achieving optimum performance.

Operating the PowerWave	36
Getting Started with Gen5	
Recommendations for Optimum Performance	37
Where to Go Next	38

Operating the PowerWave

Most users will operate the PowerWave using BioTek's Gen5 software. If you have not already done so, please follow the instructions in **Chapter 2** and the *Gen5 Getting Started Guide* for connecting the host computer and installing Gen5.

For custom applications, BioTek provides a serial communication protocol (PN 7266201-SP). Contact your BioTek dealer.

Getting Started with Gen5

These instructions briefly describe how to create and run an Experiment in Gen5. For more information, or if the instructions below do not match what you see in Gen5, refer to the **Gen5 Getting Started Guide** and help system.

For Gen5 version 2.x:

- Start Gen5.
- If the Task Manager appears, select **Read Now > New** and skip to step 4. Otherwise, select **File > New Task** from the main view.
- Select **Read Now > New**. Gen5 will open the procedure dialog. Skip to step 4.

For Gen5 version 1.x:

- Start Gen5. If the Welcome screen appears, select **Read a Plate** and skip to step 4. Otherwise, select **File > New Experiment** from the main view.
- Select **Default Protocol** and click **OK**. Gen5 will open the Experiment workspace, which includes the Protocol menu tree and Plate screen.
- Select **Plate** > **Read** or click the Read Plate icon. The Procedure dialog will open. Go to step 4.

For any version:

- Select a **Plate Type**.
- Click **Read** to open the Read Step dialog.
- 6. Select a **Read Type**.
- Select or enter the wavelength(s) at which the plate will be read.
- Define other reading parameters as desired. Click the **Help** button for assistance.
- When complete, click **OK** to return to the Procedure dialog.
- 10. Click **OK** to save and close the Procedure dialog.

- Gen5 version 1.x only: The Plate Reading dialog will open. Enter any desired information, place the plate on the carrier, and then click **READ** to begin the plate read. If the Save As dialog opens, enter a File name, choose a file location (Save in:) and click Save.
- 11. Click **OK** when the Load Plate dialog appears. The plate will be read.
 - To view the raw data results, use the Data drop-down arrow in the Plate screen to select one wavelength. The results will be displayed for the selected wavelength. Repeat, for other wavelengths.
 - To analyze, manipulate, or print results, Protocol parameters should be defined. Refer to the Gen5 Help system for instructions.

Recommendations for Optimum Performance

- Microplates should be perfectly clean and free of dust or bottom scratches. Use new microplates from sealed packages. Do not allow dust to settle on the surface of the solution; use microplate covers when not reading the plate. Filter solutions to remove particulates that could cause erroneous readings.
- Although the PowerWave supports most flat, U-bottom, and V bottom microplates, the reader achieves optimum performance with optically clear, flat-bottomed wells.
- Non uniformity in the optical density of the well bottoms can cause loss of accuracy, especially with U- and V-bottom polyvinyl microplates. Check for this by reading an empty microplate. Dual wavelength readings can eliminate this problem, or bring the variation in density readings to within acceptable limits for most measurements.
- Inaccuracy in pipetting has a large effect on measurements, especially if smaller volumes of liquid are used. For best results, use at least 100 µL per well in a 96-well plate and 25 µL in a 384-well plate (if supported).
- Dispensing solution into 384-well plates often traps air bubbles in the wells, which may result in inaccurate readings. A dual-wavelength reading method usually eliminates these inaccuracies; however, for best results, remove the air bubbles by degassing the plate in a vacuum chamber before reading.
- The inclination of the meniscus can cause loss of accuracy in some solutions, especially with small volumes. Agitate the microplate before reading to help bring this problem within acceptable limits. Use Tween® 20, if possible (or some other wetting agent) to normalize the meniscus. Some solutions develop menisci over a period of several minutes. This effect varies with the brand of microplate and the solution composition. As the center of the meniscus drops and shortens the light path, the density readings change. The meniscus shape will stabilize over time.

Where to Go Next

- Refer to the *Gen5 Getting Started Guide* and help system to learn more about using Gen5 to operate the PowerWave.
- Review the remaining chapters in this Operator's Manual to learn how to test the performance of the PowerWave, clean and maintain the reader, and troubleshoot problems.

Instrument Qualification

This chapter discusses procedures for qualifying the reader's initial and ongoing performance.

Recommended Qualification Schedule 4 Qualification Procedures 4 System Test 4 Absorbance Plate Test 4	11 11 14
Liquid Testing4	1 7

Overview

This chapter contains BioTek Instruments' recommended Installation Qualification (IQ), Operational Qualification (OQ), and Performance Qualification (PQ) procedures for the PowerWave Microplate Spectrophotometer.

Every PowerWave reader is fully tested at BioTek prior to shipment and should operate properly upon initial setup. If you suspect that a problem occurred during shipment, if you have received the equipment after returning it to the factory for service, and/or if regulatory requirements dictate that you qualify the equipment on a routine basis, you should perform the procedures outlined in this chapter.

Recommended Qualification Schedule

This schedule defines the factory-recommended intervals for qualifying a PowerWave used two to five days a week. The actual frequency, however, may be adjusted depending on your usage of the instrument. This schedule assumes the reader is properly maintained as outlined in *Chapter 5*.

❖ The risk factors associated with your tests may require that the Operational and Performance Qualification procedures be performed more or less frequently than indicated here.

	Installation Qualification	Operational Qualification		rmance fication
	Initially	Initially/ Annually	Monthly	Quarterly
System Test	✓	✓	✓	
Absorbance Plate Test		✓	✓	
Liquid Test 1*		✓		✓
Liquid Test 2*		✓		✓
Liquid Test 3**		✓		✓

Run Liquid Test 1 if you have an Absorbance Test Plate. If you do not have a Test Plate, run Liquid Test 2.

^{**} Liquid Test 3 is optional; it is provided for sites requiring verification at wavelengths lower than those attainable with the Absorbance Test Plate.

Qualification Procedures

Your reader was fully tested at BioTek prior to shipment and should operate properly upon initial setup. If you suspect that problems occurred during shipment or if regulatory requirements dictate that Operational and/or Performance Qualification is necessary, you should perform the following tests.

- **System Test:** Verifies proper gains, bulb operation, low electronic noise, and incubator functionality. The test report includes the reader's serial number and basecode software part number and version number.
- **Absorbance Plate Test:** Uses BioTek's Absorbance Test Plate to confirm the mechanical alignment, optical accuracy/linearity, repeatability, channel-to-channel variation, and wavelength accuracy of the instrument.
- **Liquid Tests:** Uses liquid solutions in a microplate to confirm mechanical alignment, optical accuracy/linearity, repeatability, and channel-to-channel variation of the instrument.

System Test

The PowerWave automatically runs the System Test each time it is powered on. The test can also be run (and a report generated) using Gen5. The System Test report should be printed to document periodic testing and for troubleshooting purposes. See a sample System Test report on page 42.

If the power-up System Test fails, the reader will beep. If this happens, press the carrier eject button to stop the beep and then attempt to run the System Test using Gen5 to retrieve an error code from the instrument.

See *Chapter 6* for a list of possible error codes and their probable causes.

- 1. Adjust the wavelength table, if desired:
 - From Gen5's main screen, select **System > Reader Configuration**.
 - Highlight the PowerWave and click **View/Modify > Setup**.
 - Click the Absorbance tab.
 - Enter the desired wavelength values and click **Send Wavelengths** to download them to the reader.
- 2. Run the test:
 - Select System > Diagnostics > Run System Test.
 - The test will run and results will appear in a pass/fail format.

Sample System Test Report

Test Date/Time: 11/06/2009 10:33:27 AM

Reader: 123456
Operator: System Administrator
Comments: System Test performed during IQ

Operator ID:__

Notes:__

			SYST	EM SELF	TEST				
7280201 Ve	rsion 1	.20	1287	90					
Lambda: 200		3.28	Resets						
Channel:	Ref	1	2	3	4	_	6	7	8
Air:	20506	35689	37068	35775	39588	36708	32962	35604	34634
Dark:	9911	10002	10031		9980	9904	9930	9869	9974
Delta:	10595	25687	27037	25895	29608	26804	23032	25735	24660
Lambda: 405	Gain:	2.53	Resets	: 2					
Channel:	Ref	1	2	3	4	5	6	7	8
Air:	20072	35592	39355	36205	38232	37636	33564	36682	36095
Dark:	9914	9960	9975	9900	9945	9915	9928	9895	9948
Delta:	10158	25632	29380	26305	28287	27721	23636	26787	26147
Lambda: 490	Gain:	1.72	Resets	: 2					
Channel:	Ref	1	2	3	4	5	6	7	8
Air:	20280	35642	39779			37962	33494	36949	36219
Dark:	9919	9949	9960	9909	9939	9919	9928	9905	9942
	10361	25693	29819	26401	28547	28043	23566	27044	26277
Lambda: 550	Gain:	2.51	Resets	: 1					
Channel:	Ref	1	2	3	4	5	6	7	8
Air:	20248	35530	39702		38344	37862	33339	36781	36154
Dark:	9914	9985	10007	9891	9967	9910	9930	9881	9963
Delta:	10334	25545	29695	26208	28377	27952	23409	26900	26191
Lambda: 630	Gain:	3.94	Resets	: 1					
Channel:	Ref	1	2	3	4	5	6	7	8
Air:	20230	35302	39381	35837	38113	37659	33193	36503	35852
Dark:	9907	10018	10053	9870	9990	9900	9931	9854	9984
Delta:	10323	25284	29328	25967	28123	27759	23262	26649	25868
Lambda: 999	Gain:	5.33	Resets	: 1					
Channel:	Ref	1	2	3	4	5	6	7	8
Air:	19977	34186	39376	35131	37500	36779	32458	36029	35249
Dark:	9900	10050	10095	9849	10010	9888	9934	9831	10004
Delta:	10077	24136	29281	25282	27490	26891	22524	26198	25245

	Ref 9878 9877 1			3 9847 9846 1			6 9919 9917 2	7 9836 9835 1		
INCUBATOR SEI	LF TEST									
Temperature S	Setpoin	t: 37.	0	Current	Average	: 37.1	A	/D Test	: PASS	
Zone 1: 37.1 Zone 2: 37.0 Zone 3: 37.1 Zone 4: 37.1	Min: Min:	37.0 37.1	Max: Max:		Range:	PASS PASS PASS PASS	Therr Therr	mistor: mistor: mistor: mistor:	PASS PASS	
Upper Left Co Lower Left Co Lower Right (Upper Right (Delta 1: 167	orner: Corner: Corner:	x= 24 x= 111 x= 103	48 y: 40 y:	= 9220 = 8944 = 8912 = 9172	the F	PASS or	TOR SELF FAIL con- zones in	dition of	the	
Delta 2: 1036 Delta 3: 917 Delta 4: 891	72 - 92	20= -		last au instrur	CAL ANAL' utocalibrat ment. AU1	ion proc TOCAL is	edure pe not inclu	rformed uded in th	on the	
Middle Sensor Tested: Delta:	c:	x= 138 138		SYSTE	M TEST re	esuit giv	en below			
SYSTEM TEST I	PASS 🔨			M TEST resto the opt				FAIL and		

Absorbance Plate Test

This test uses BioTek's Absorbance Test Plate to confirm the mechanical alignment, optical density accuracy/linearity, repeatability, channel-to-channel variation and wavelength accuracy of the PowerWave. The test compares the reader's optical density and wavelength measurements to NIST-traceable values.

❖ If you do not have an absorbance test plate, you can run Liquid Test 2 to test accuracy/linearity, repeatability, and alignment. See page 50.

Test Plate and Certificates

To run the Absorbance Plate Test, you will need BioTek's Absorbance Test Plate (PN 7260522), with its accompanying certificates.

- The **Standards Certificate** contains standard OD values for the filters at several different wavelengths.
- The **Peak Wavelength Certificate** contains one or more peak wavelength values for the glass filter in position C6 on the plate. Each value has a valid test range associated with it. For example, a Peak Wavelength value may be 586 nm with a test range of 580 to 590 nm (or tolerance values of -6/+4).
- The instructions provided below and on the following page are guidelines. Refer to the Gen5 Help system for more information.

Enter the Absorbance Test Plate Data

Before performing the Absorbance Plate Test for the first time (and after the plate is recalibrated by BioTek), enter data from the Standards and Peak Wavelengths Certificates into Gen5:

Select System > Diagnostics > Test Plates > Add/Modify Plates and click **Add**. Click **Help** for guidance when setting the wavelengths and entering the OD and peak wavelength value(s).

Run the Absorbance Plate Test

Place the Absorbance Test Plate in the carrier so that well A1 is in the right rear corner of the carrier.

- Select System > Diagnostics > Test Plates > Run.
- If the Instrument Selection dialog appears, select the appropriate reader.
- If the Select Test Plate dialog appears, select the appropriate test plate.
- Enable 'Perform peak wavelength test' and enter any required information
- Click Start Test.

Sample Absorbance Plate Test Report

```
Universal Test Plate Analysis
Wavelength: 405
Alignment Results
A1= 0.000 PASS
             A12= 0.000 PASS H1= 0.000 PASS H12= 0.000 PASS
Accuracy Results
             C01
                   D04
                         E02
                               F05
                                      G03
                                            H06
             0.147 2.945 0.618 2.279 1.133 1.701
  STANDARD
            0.139 2.914 0.613 2.265 1.128 1.694
  DATA
  RESULT
             PASS PASS PASS PASS
                                           PASS
             F12
                   E09 D11
                               C08
                                     B10
                                           A07
  STANDARD 0.147 2.945 0.618 2.279 1.133 1.701
             0.145 2.910 0.613 2.264 1.125 1.692
  DATA
  RESULT
            PASS
                  PASS
                        PASS
                               PASS
                                     PASS
                                           PASS
Repeatability Results
             C01
                   D04
                         E02
                               F05
                                     G03
                                            Н06
  READ1
             0.139 2.914 0.613 2.265 1.128 1.694
             0.141 2.913 0.616 2.265 1.128 1.694
  READ2
  RESULT
             PASS
                  PASS PASS
                               PASS PASS
                                           PASS
             F12
                   E09
                         D11
                               C08
                                     B10
                                            A07
  READ1
            0.145 2.910 0.613 2.264 1.125 1.692
           0.140 2.914 0.614 2.265 1.126 1.692
  READ2
  RESULT
            PASS PASS PASS PASS PASS
Turnaround Results
             C01
                   D04
                        E02
                               F05
                                     G03
                                            H06
             0.139 2.914 0.613 2.265 1.128 1.694
  READ1
             F12
                   E09 D11
                               C08
                                     B10
                                           A07
             0.145 2.910 0.613 2.264 1.125 1.692
  READ2
  RESULT
            PASS PASS PASS PASS PASS
Spectral Scan Results
  580nm= 1.933 586nm= 2.892
  581nm= 1.956 587nm= 2.753
582nm= 2.059 588nm= 2.516
  583nm= 2.253 589nm= 2.334
  584nm= 2.511
               590nm= 2.186
  585nm= 2.785
  Test Plate Standard=587nm
  Calculated Peak=586nm PASS
```

Test Results

For the Accuracy, Linearity, Repeatability, and Channel-to-Channel Variation tests, there may not be a pass/fail indication for filter values that are beyond the specified accuracy range of the instrument.

The Absorbance Test Plate Analysis contains results for the following:

- **Mechanical Alignment**: This portion of the test measures the alignment of the microplate carrier with the optical path. A reading greater than 0.015 OD represents an out-of-alignment condition.
- **Accuracy/Linearity**: Accuracy is a measure of the absorbance (optical density) of Test Plate wells C01, D04, E02, F05, G03, and H06 as compared with known standard values contained in the plate's Standards Certificate. Linearity of the optical density readings is confirmed by default if the readings are accurate. To further verify this, perform a regression analysis on the Test Plate OD values in a program such as Microsoft® Excel. An R Square value of at least 0.99 is expected.
- **Repeatability**: Repeatability is a measure of the instrument's ability to read the same well with minimum variation between the two reads with the well in the same location.
- Channel-to-Channel Variation: This test ensures that selected channels read the same value for a filter as their paired channel when the plate is rotated 180° in the plate carrier. The channel/well "pairs" for the turnaround test are: C01/F12; D04/E09; E02/D11; F05/C08; G03/B10; H06/A07.
- Wavelength Accuracy: If 'Perform peak wavelength test' is enabled as part of the Absorbance Plate Test, the C6 filter is scanned across a specified wavelength range in 1-nm increments. The wavelength of the maximum absorbance is compared with the peak wavelength value(s) entered in the software. The accuracy of the wavelength should be ± 3 nm (± 2 nm instrument, ± 1 nm filter allowance).

For example, if the test range is 580 to 590 nm, the Certificate value is 587 nm, and the reader reports a peak value of 590 nm, then the reader meets specifications. If the reader reports 591 nm, then the reader does not meet specifications.

Liquid Testing

Conducting Liquid Tests confirms the reader's ability to perform to specification with liquid samples. Liquid testing differs from testing with the Absorbance Test Plate in that liquid in the wells has a meniscus, whereas the Test Plate's neutral density glass filters do not. The optics characteristics may differ in these two cases, thus alerting the operator to different types of problems.

- If you have the Absorbance Test Plate, you only need to perform **Liquid Test 1** for routine testing.
- If you do not have the Absorbance Test Plate, you can test the accuracy/linearity, repeatability, and alignment of the reader by performing Liquid Test 2.
- If you must test the reader's performance at 340 nm, perform **Liquid Test 3**.
- BioTek offers a dye solution (PN 7120779, 25 mL; or 7120782, 125 mL) that may be used in the stock solution formulation for Liquid Tests 1 and 2, or, if you prefer, you may use Solution A described on the next page. The purpose of the formulation is to create a solution that absorbs light at approximately 2.000 OD full strength when dispensed at 200 µL in a flat-bottom microplate well.
- Alternatively, any solution that gives a stable color will suffice. (This includes substrates incubated with an enzyme preparation and then stopped with an acidic or basic solution.) Some enzyme/substrate combinations that may be used as alternates to the described dye are shown below.

Typical Enzy	me-Substrate	Combinations	and Stopping	1 Solutions
I ypicai Liiz	ymic-Jubstrate	Combinations		Jointions

Enzyme	Substrate	Stopping Solution
Alkaline Phosphate	o-nitrophenyl phosphate	3N sodium hydroxide
beta-Galactosidase	o-nitrophenyl -beta-D galactopyranoside	1M sodium carbonate
Peroxidase	2,2'-Azino di-ethylbenzothiazoline- sulfonic acid (ABTS)	citrate-phosphate buffer, pH 2.8
Peroxidase	o-phenylenediamine	0.03N sulfuric acid

Absorbance Liquid Test 1

Absorbance Liquid Test 1 confirms repeatability and alignment of the reader when a solution is used in the microplate. If these tests pass, then the lens placement and optical system cleanliness are proven.

Materials

Manufacturer part numbers are subject to change.

- New 96-well, clear, flat-bottom microplate (Corning Costar® #3590 recommended)
- Stock Solution A or B, which may be formulated by diluting a dye solution available from BioTek (A) or from the ingredients listed below (B).

Solution A

- BioTek QC Check Solution No. 1 (PN 7120779, 25 mL; PN 7120782, 125 mL)
- Deionized water
- 5-mL Class A volumetric pipette
- 100-mL volumetric flask
- Pipette a 5-mL aliquot of BioTek QC Check Solution No. 1 into a 100-mL volumetric flask.
- Add 95 mL of DI water; cap and shake well. The solution should measure approximately 2.000 OD when using 200 µL in a flat-bottom microwell.

Solution B

- Deionized water
- FD&C Yellow No. 5 dye powder (typically 90% pure)
- Tween® 20 (polyoxyethylene (20) sorbitan monolaurate) or BioTek wetting agent (PN 7773002) (a 10% Tween® solution)
- Precision balance with capacity of 100g minimum and readability of 0.001g
- Weigh boat
- 1-liter volumetric flask
- Weigh out 0.092 g of FD&C Yellow No. 5 dye powder into a weigh boat.
- Rinse the contents into a 1-liter volumetric flask. 2.
- 3. Add 0.5 mL of Tween 20, or 5 mL of BioTek's wetting agent.
- Fill to 1 liter with DI water; cap and shake well. The solution should measure approximately 2.000 OD when using 200 µL in a flat-bottom microwell.

Prepare the Plate

- ❖ A 96-well, flat-bottom microplate is required for this test (Corning Costar® #3590 is recommended). Use a new microplate; any fingerprints or scratches may cause variations in readings.
 - 1. Using a freshly prepared stock solution (Solution A or B), prepare a 1:2 dilution using deionized water (one part stock, one part deionized water; the resulting solution is a 1:2 dilution).
 - Pipette 200 μL of the concentrated solution into the first column of wells in the new 96-well microplate.
 - Pipette 200 µL of the diluted solution into the second column of wells.



Important! After pipetting the diluted test solution into the microplate and *before* reading the plate, we strongly recommend shaking the plate at Variable speed for four minutes. This will allow any air bubbles in the solution to settle and the meniscus to stabilize. Alternatively, wait 20 minutes after pipetting the diluted test solution before reading the plate.

Read the Plate

- Read the microplate five times at 405 nm using normal reading mode, single wavelength, no blanking ("Normal" plate position).
- Rotate the microplate 180° so that well A1 is now in the H12 position. Read the plate five more times ("Turnaround" plate position).
- Print out the ten sets of raw data, or export them to an Excel spreadsheet.

Calculations

- Calculate the mean value for each physical well location in columns 1 and 2 for the five plates read in the Normal position, and then again for the five plates read in the Turnaround position. This will result in 32 mean values.
- 2. Perform a mathematical comparison of the mean values for each microwell in its Normal and Turnaround positions (A1/H12, A2/H11, B1/G12, B2/G11, and so on). In order to pass this test, the differences in the compared mean values must be within the accuracy specification for the instrument.

For Example:

If the mean value for well A1 in the Normal position is 1.902, where the specified accuracy is $\pm 1.0\% \pm 0.010$ OD, then the expected range for the

mean of the same well in its Turnaround (H12) position is 1.873 to 1.931 OD.

```
1.902 * 0.01 + 0.010 = 0.029; 1.902 - 0.029 = 1.873; 1.902 + 0.029 = 1.931
```

If any set of mean values is out of the expected range, review the other three sets of mean values for the same channel pair. For example, if the A1/H12 comparison fails (the wells are not within the expected range of each other), review the comparisons of A2/H11, H1/A12, and H2/A11. If two or more sets of mean values for a channel pair are out of the expected range, there is a problem with one of the eight read channels. If only one of the four mean values results in a failure, check the well for debris and the plate for scratches or fingerprints.

Accuracy Specification:

For comparison in this test, the following accuracy specifications are applied, using Normal/Standard read mode and a 96-well microplate.

```
\pm 1.0\% \pm 0.010 OD from 0.000 to 2.000 OD
\pm 3.0\% \pm 0.010 OD from 2.000 OD to 3.000 OD
```

Absorbance Liquid Test 2

The recommended method of testing the instrument performance is to use the Absorbance Test Plate to confirm alignment, accuracy/linearity, and repeatability. If a Test Plate is not available, Liquid Test 2 can be used for these tests.

Materials

- A new 96-well, flat bottom microplate (Corning Costar® #3590 is recommended)
- Ten test tubes, numbered consecutively, stored in a rack
- Calibrated hand pipette (Class A volumetric pipette recommended)
- Solution A or B (these are the same solutions as for Liquid Test 1)
- A 0.05% solution of deionized water and Tween® 20

Prepare the Dilutions

Create a percentage of dilution series, beginning with 100% of the original concentrated stock solution (A or B) in the first tube, 90% of the original solution in the second tube, 80% in the third tube, all the way to 10% in the tenth tube. Dilute using the 0.05% solution of deionized water and Tween® 20. This solution can also be made by diluting the BioTek wetting agent 200:1.

Tube Number	1	2	3	4	5	6	7	8	9	10
Volume of Original Solution (mL)	20	18	16	14	12	10	8	6	4	2
Volume of 0.05% Tween solution (mL)	0	2	4	6	8	10	12	14	16	18
Absorbance expected if original solution is 2.000 OD at 200 µL	2.0	1.8	1.6	1.4	1.2	1.0	0.8	0.6	0.4	0.2

Test Tube Dilutions

❖ The choice of dilutions and the absorbance of the original solution can be varied. Use this table as a model for calculating the expected absorbances of a series of dilutions, given a different absorbance of the original solution.

Prepare the Plate

- 1. Pipette $200 \mu L$ of the concentrated solution from tube 1 into each well of the first column, A1 to H1, of a new flat-bottom microplate.
- 2. Pipette 200 μL from each of the remaining tubes into the wells of the corresponding column of the microplate (tube 2 into wells A2 to H2, etc.).
 - After pipetting the diluted test solution into the microplate and before reading the plate, we strongly recommend shaking the plate at Variable speed for four minutes. This will allow any air bubbles in the solution to settle and the meniscus to stabilize. Alternatively, wait 20 minutes after pipetting the diluted test solution before reading the plate.

Linearity and Repeatability Tests

- 1. Read the microplate prepared above using a normal mode dual wavelength at 450 nm with 630 nm as the blank. Repeat the read four times for a total of five reads. Retain the plate for the Alignment test.
- 2. Print out the five sets of Delta OD data, or export them to an Excel spreadsheet.
- 3. Calculate the results for Linearity:
 - Calculate the mean absorbance for each well, and average the means for each concentration.
 - Perform a regression analysis on the data to determine if there is adequate linearity.

Since it is somewhat difficult to achieve high pipetting accuracy when conducting linear dilutions, an R Square value of at least 0.99 is considered adequate.

- 4. Calculate the results for Repeatability:
 - Calculate the mean and standard deviation for the five readings taken in Step 1 at each concentration. Only one row of data needs to be analyzed.
 - For each mean below 2.000 OD, calculate the allowed deviation using the repeatability specification for a 96-well plate of $\pm 1.0\% \pm 0.005$ OD. If above 2.000 OD, apply the $\pm 3.0\% \pm 0.005$ specification.
 - The standard deviation for each set of readings should be less than the allowed deviation.

Example: Absorbance readings of 1.950, 1.948, 1.955, 1.952, and 1.950 will result in a mean of 1.951, and a standard deviation of 0.0026. The mean (1.951) multiplied by 1% (1.951 * 0.010) = 0.0195, which, when added to the 0.005 (0.0195 + 0.005) = 0.0245 OD, which is the allowable deviation. Since the standard deviation is less than this value, the reader meets the test criteria.

Repeatability Specification:

 $\pm 1\% \pm 0.005$ OD from 0.000 OD to 2.000 OD

 $\pm 3\% \pm 0.005$ OD from 2.000 OD to 3.000 OD

Channel-to-Channel Variation and Alignment

- Using the plate prepared for the tests above, conduct a turnaround test by reading the plate with the A1 well in the H12 position five times. This test results in four comparisons of each channel to its corresponding channel, two in column 1, and two in column 2.
- Calculate the means of the wells in columns 1 and 2 in the normal plate position (data is from the tests above) and in the turnaround position (from Step 1 above). Compare the mean reading for well A1 to its mean reading when in the H12 position.
- Compare the mean values for the other wells to their corresponding mean values with the well in the turnaround position. (Compare B1 to G12, C1 to F12, D1 to E12, E1 to D12, F1 to C12, G1 to B12, H1 to A12, A2 to H11, and B2 to G11, etc.). The difference in the values for any two corresponding wells should be within the accuracy specification for the instrument.

For example: If the mean of well A1 in the normal position is 1.902, where the specified accuracy is $\pm 1.0\% \pm 0.010$ OD, then the expected range for the mean of the same well in the H12 position is 1.873 to 1.931 OD. (1.902 * 1% = 1.873)0.019 + 0.010 = 0.029, which is added and subtracted from 1.902 for the range.)

If any set of well values is out of the expected range, review the other three sets for the same channel pair. Thus, if A1 and H12 are not within range of each other, review the compliance of H1 to A12, A2 to H11, and H2 to A11.

This will confirm that there is a problem in one of the eight read channels, or indicate that the result of one set of wells was in error. If any two sets of well values for a channel pair are out of the allowed accuracy range, there may be contamination on, or a problem with, one of the lenses.

If the four corner wells are within the accuracy range, the reader is also in alignment.

Absorbance Liquid Test 3 (Optional)

Materials

- A new 96-well, flat-bottom microplate (Corning Costar® #3590 is recommended)
- Calibrated hand pipette(s)
- Beakers and graduated cylinder
- Precision balance with a readability of 0.01 g
- Buffer solution as described below

Buffer Solution

- Deionized water
- Phosphate-buffered saline (PBS), pH 7.2-7.6, Sigma® tablets #P4417 (or equivalent)
- β-NADH Powder (β-Nicotinamide Adenine Dinucleotide, Reduced Form) Sigma bulk catalog number N 8129, or preweighed 10-mg vials, Sigma number N6785-10VL (or BioTek PN 98233). Store the powder according to the guidelines on its packaging.
- 1. Prepare a PBS solution using Sigma® tablets.
- In a beaker, mix 50 mL of the PBS solution with 10 mg of the β -NADH powder and mix thoroughly. This is the **100% Test Solution**.

Prepare the Plate

- Prepare the 75% Test Solution by mixing 15 mL of the 100% Test Solution with 5 mL of the PBS solution.
- Prepare the 50% Test Solution by mixing 10 mL of the 100% Test Solution with 10 mL of the PBS solution.
- Pipette the three solutions into the new 96-well microplate: 3.
 - > 150 μL of the 100% Test Solution into all wells of columns 1 and 2
 - \blacktriangleright 150 μ L of the 75% Test Solution into all wells of columns 3 and 4
 - > 150 μL of the 50% Test Solution into all wells of columns 5 and 6



Important! After pipetting the diluted test solution into the microplate and *before* reading the plate, we strongly recommend shaking the plate at Variable speed for four minutes. This will allow any air bubbles in the solution to settle and the meniscus to stabilize. Alternatively, wait 20 minutes after pipetting the diluted test solution before reading the plate.

Read the Plate

- Read the microplate five times using Normal mode, single wavelength at 340 nm, no blanking (or blank on air).
- Print the five sets of raw data or export it to an Excel spreadsheet using Gen5.

Analyze the Results

- 1. For each well, calculate the mean and standard deviation of the five readings.
- For each mean calculated in step 1, calculate the allowed deviation using the repeatability specification for a 96-well plate of ± 1.0% ± 0.005 OD (mean * 0.01 + 0.005).
- For each well, compare the standard deviation calculated in step 1 with the allowed deviation calculated in step 2. The standard deviation should be less than the allowed deviation.

Example: Five readings in well A1 of 0.802, 0.802, 0.799, 0.798, and 0.801 will result in a mean of 0.8004, and a standard deviation of 0.0018. The mean multiplied by 1.0% (0.8004 * 0.010) = 0.008, which, when added to the 0.005 (0.008) +0.005) = 0.013, which is the allowable deviation for well A1. Since the standard deviation for well A1 is less than 0.013, the reader meets the test criteria.

- Calculate the results for Linearity:
 - For each of the three dye concentrations, calculate the mean absorbance for the wells containing that solution (mean of wells A1 to H2, A3 to H4, and A5 to H6).
 - Perform a regression analysis on the data to determine if there is adequate linearity.

Expected Results: Since it is somewhat difficult to achieve high pipetting accuracy when conducting linear dilutions, an R-Square value greater than or equal to 0.990 is considered adequate.

Chapter 5 Decontamination

This chapter contains procedures for decontaminating and cleaning the PowerWave.

Purpose	58
Clean Plate Carrier and Exposed Surfaces	58
Decontamination	59

Purpose

Any laboratory instrument that has been used for research or clinical analysis is considered a biohazard and requires decontamination prior to handling. Intact skin is generally considered an effective barrier against infectious organisms; however, small abrasions and cuts may not always be visible.

Decontamination minimizes the risk to all who come in contact with the instrument during shipping, handling, and servicing. Decontamination is required by the U.S. Department of Transportation regulations.

Persons performing the decontamination process must be familiar with the basic setup and operation of the instrument.



BioTek Instruments, Inc. recommends the use of the following decontamination solutions and methods based on our knowledge of the instrument and recommendations of the Centers for Disease Control and Prevention (CDC). Neither BioTek nor the CDC assumes any liability for the adequacy of these solutions and methods. Each laboratory must ensure that decontamination procedures are adequate for the Biohazard(s) they handle.



Warning! Internal Voltage. Turn off and disconnect the PowerWave from its power supply for all cleaning and decontamination operations.



Warning! Wear prophylactic gloves when handling contaminated instruments. Gloved hands should be considered contaminated at all times; keep gloved hands away from eyes, mouth, nose, and ears. Eating and drinking while decontaminating instruments is not advised.



Warning! Mucous membranes are considered prime entry routes for infectious agents. Wear eye protection and a surgical mask when there is a possibility of aerosol contamination. Intact skin is generally considered an effective barrier against infectious organisms; however, small abrasions and cuts may not always be visible. Wear protective gloves when handling contaminated instruments.



Important! Do not immerse the instrument, spray it with liquid, or use a "wet" cloth. Do not allow water or other cleaning solution to run into the interior of the instrument. If this happens, contact BioTek's Technical Assistance Center.

Clean Plate Carrier and Exposed Surfaces



Important! Turn off the PowerWave and disconnect it from the power supply for the cleaning procedure.

A regular cleaning regimen is recommended to keep the instrument free of dust and particulates that can cause erroneous readings. Exposed surfaces may be cleaned (not decontaminated) with a cloth moistened (not soaked) with water or water and a mild detergent. You will need:

- Deionized or distilled water
- Clean lint-free cotton cloths
- Mild detergent (optional)
- Turn on the PowerWave and press the carrier eject button to eject the microplate carrier.
- Turn off and unplug the reader from the power supply.
- Moisten a clean, lint-free cloth with water, or with water and the mild detergent. **Do** not soak the cloth.
- Wipe the plate carrier and all exposed surfaces of the instrument.
- 5. If detergent was used, wipe all surfaces with a cloth moistened with water.
- Use a clean, dry lint-free cloth to dry all wet surfaces.

Decontamination

Tools and Supplies

- Sodium hypochlorite (NaClO, or bleach)
- Deionized or distilled water
- Safety glasses
- Surgical mask
- Protective gloves
- Lab coat
- Biohazard trash bags
- 125-mL beakers
- Clean, lint-free cotton cloths



Warning! The bleach solution is caustic; wear gloves and eye protection when handling the solution.

- Turn on the PowerWave and press the carrier eject button to eject the carrier. 1.
- 2. Turn off and unplug the instrument from the power supply.
- 3. Prepare an aqueous solution of 0.5% sodium hypochlorite (NaClO, or bleach).
 - Check the % NaClO of the bleach you are using; this information is printed on the side of the bottle. Commercial bleach is typically 10% NaClO; if this is the case, prepare a 1:20 dilution. Household bleach is typically 5% NaClO; if this is the case, prepare a 1:10 dilution.
- Moisten a clean, lint-free cloth with the bleach solution. Do not soak the cloth. 4.
- 5. Wipe the plate carrier and all exposed surfaces of the instrument.
- Allow the instrument to dry for 20 minutes for thorough decontamination by the bleach.
- 7. Moisten a cloth with deionized or distilled water and wipe all surfaces of the instrument that have been cleaned with the bleach solution.
- 8. Use a clean, dry lint-free cloth to dry all wet surfaces.
- Discard the used gloves and cloths, using a Biohazard trash bag and an approved Biohazard container.

Chapter 6 Troubleshooting

This chapter lists error codes that may appear during operation of the PowerWave, and provides troubleshooting tips.

Error Codes	62
General Errors	63
Fatal Errors (A###)	67

Error Codes

An error code is displayed in the software as a four-digit identifier. The first digit will be 0, 1, 2, 3, or A.

- 0, 1, 2, or 3 denote a noncritical error, which means it is still possible for the PowerWave to communicate with the controlling software and run a reader system test (see below). See General Errors starting on page 63.
- "A" denotes a more serious error with the memory or processing, which requires the reader to be turned off/on before any diagnostics can be performed. If the error reappears, contact BioTek TAC for troubleshooting assistance (see page 23). See Fatal Errors on page 67.

If an error code is displayed in Gen5, run a System Self-Test for diagnostic purposes:

• From the main screen, select **System > Diagnostics > Run System Test**.

Error Codes During Operation with the BioStack

Error codes may appear in Gen5 during operation of the PowerWave with the BioStack Microplate Stacker.

Gen5 error codes display in a negative value format, for example: -8, -100.

Refer to the Bio Stack Operator's Manual for a list of error codes and their descriptions.

See Product Support and Service in Chapter 1 for contact information for BioTek's Technical Assistance Center ("TAC").

General Errors

Code	Description and Probable Causes
0200	X-axis opto sensor failed to transition.
	This error indicates that a motor was not able to move to its "home" position, as registered by feedback from an optical sensor, or it failed to transition after moving away from the home position.
	Dirty x-axis rail or dry bearings are causing too much friction.
	Defective or broken optical sensor.
	Defective motor controller PCB.
0201	Order sorting (bandpass) filter wheel did not home. • Filter wheel is loose.
	Filter wheel is obstructed by too close proximity to the motor gear.
	Defective or broken optical sensor.
	Defective motor, motor controller PCB, or cable.
0202	Y-axis opto sensor failed to transition.
	This error indicates that a motor was not able to move to its "home" position, as registered by feedback from an optical sensor, or it failed to transition after moving away from the home position.
	Dirty y-axis rail or dry bearings are causing too much friction.
	Defective or broken optical sensor.
	Defective motor controller PCB.
	Note: In cases where a sensor is not functioning, the motor will drive the axis to its mechanical stop and generate substantial noise.
0303	Monochromator did not find home.
	During the instrument initialization, the monochromator is homed by rotating the monochromator mirror until the white light (full light) is detected. This requires a fully functional flash lamp/detection system. • Defective analog PCB.
	Defective flash lamp and or flash lamp power supply (inconsistent flashes) (high probability).
	Defective motor/power PCB.
	Defective monochromator (low probability).
0400	Carrier x-axis failed positional verify.
	Motor x-axis failed to reach the same position when moved a known number of steps from the home position and back.
	Dirty rail or dry bearings are causing too much friction.
0401	Order sorting (bandpass) filter wheel failed positional verify.
	Filter wheel obstructed by motor gear.
	Motor gear loose on motor shaft.
0402	Carrier y-axis failed positional verify.
	Dirty rail or dry bearings are causing too much friction.

Code	Description and Probable Causes
0403	Monochromator failed to find the zero order position (white light).
	The order sorting (bandpass) filter wheel is homed and moved to the open hole position. The monochromator is moved until the optical system detects saturation (home). It is ther moved to a known number of steps away from home and then moved back the same number of steps, expecting to see light saturation point. The error is indicating the saturation did not clear or appear.
	Flash lamp is missing flashes or is not flashing.
	The optic system does not detect the saturation.
	Defective monochromator.
0500	Measurement or reference channel is saturated during a spectral scan.
	This error indicates the light signal level in one of the channels reached 65,535 counts during Lambda calibration within the spectral scan.
0501- 0508	Measurement or reference channel is saturated during a spectral scan. This error indicates the light signal level in the channel indicated by the last digit (1 through 8) of the error code reached 65,535 counts during the spectral scan or calibration The lamp is not properly aligned and there is too much light. The A/D reference voltage is not at the 4.5 V. The analog PCB is defective.
	Order sorting filter has degraded.
0503	Monochromator failed positional verify due to saturation.
0303	This error indicates that, during initialization, the monochromator failed positional verify, channel 3 failed during calibration or spectral scan.
0511-	Measurement channel A/D signal saturated.
0568	This error indicates the light signal level reached 65,535 counts for one of the lambda values in the table after calibration, prior to a read, or during a read or optics test.
0600	Gain out of range for the target air readings. Reference channel = hot channel.
	This error indicates that the measurement channel signal gain is out of range necessary to ensure the reader's performance to specifications. During reader calibration, the gain selected is 36.56. • Flash lamp • Monochromator
	Lamp power supply
0701-	Channel failed noise test greater than 20 counts during optics test.
0708	This error indicates significant variations in background electronic noise were detected, when blocking the light and increasing the gain to maximum.
	• Electrical noise may be penetrating the measurement chamber. The bottom and top shrouds are part of the electrical shielding.
	The coaxial cable ground may be floating or disconnected.
	There may be an ambient light leak. Ensure that the plate carrier door is properly closed Applied BCB failure, poissy photo detector.
	 Analog PCB failure; noisy photo-detector. Internal electronic noise may be caused by a faulty analog PCB or faulty internal
	grounding.
0801-	Channel failed noise offset < 10 and > 2000 during optics test.
8080	This error indicates that background electronic signal detected is outside of acceptable limits at maximum gain when blocking the light.
	The photo-detector is not connected or is defective, yielding a noise reading of zero.
	The photo-detector is too noisy and is defective.
0901- 0908	Channel dark range is < 100 or > 20000 during calibration, or < 100 during a filter test.

Code	Description and Probable Causes
0911- 0968	Measurement channel dark range is < 100 during a read in enhanced mode (64 flashes), or prior to a read or optics test. The reference channel dark current value has changed since the last optics test measurement by more than 10%, or the dark value is less than 100. The last number in the error code is the channel number used during the failure.
	 The photo-detector is more sensitive to temperature changes. Ambient light leakage during the read.
0A00- 0A68	Measurement channel air / blank out of range prior to a read. This error is indicating the air reading at the time of the plate read was 50% less than the air reading at the time of the optic test. The last number of the error code represents the channel at the time of failure.
0A10	Reference channel air / blank out of range for the first wavelength in a scan, or filter / reference channel air / blank out of range prior to a read. This error is indicating the air reading at the time of the plate read was 50% less than the air reading at the time of the optic test. The last number of the error code represents the channel at the time of failure. • Flash lamp has missed flashes during the read. • Dirty optics or spilled substance on the optics.
0F00- 0F08	Channel Delta out of range during calibration. The Delta of the air / dark is out of range during the calibration at a wavelength reference channel < 500 or measurement channel < 8000.
0F10	Reference channel out of range during a spectral scan. Reference channel < 500 during a spectral scan and only checking the first wavelength.
0F10- 0F60	Reference channel out of range during a read. Reference channel < 500 during a read. This error indicates that the reading has failed. The last number of the error code represents the channel.
0F10- 0F68	Channel out of range prior to a read. Reference channel < 500 or measurement channel < 8000 during an optics test or prior to a read. Reference channel out of range (< 50% or > 200%). • The flash lamp may be out of alignment. • The order sorting (bandpass) filter is degraded, and does not allow enough light energy to pass through. • Damaged reference channel optic spray. • The reference channel photodiode detection circuit is defective.
1100- 1101	Failed configuration checksum test for reader protocol or system configuration prior to a read or optics test. Last digit can be either a 0 or 1. • The flash memory on the PCB is defective or corrupt. The basecode software and/or assays may need to be re-downloaded.
1200	Lambda calibration data missing prior to a spectral scan, or autocalibration data is missing. The instrument calibration values are not loaded into the flash memory. • The PCB was changed and the flash memory does not have the calibration values loaded. • Failure in Main PCB memory. Contact BioTek TAC for more information.
1201- 1206	Lambda table calibration data missing from reader. This error occurs when the controlling PC requests the Lambda wavelength calibration data, and one of the wavelengths does not have calibration data in memory (not calibrated). The last digit represents the Lambda.

Code	Description and Probable Causes
1300	Carrier not homed in the x-axis.
	This error is only seen if an error 0200 is ignored. See the probable causes for 0200.
1301	Order sorting (bandpass) filter wheel not homed.
	This error is only seen if an error 0201 is ignored. See the probable causes for 0201.
1302	Carrier not homed in the y-axis.
	This error is only seen if an error 0202 is ignored. See the probable causes for 0202.
1501- 1504	Temperature zone out of range (the last digit is the zone number failing).
1511- 1514	Thermistor failed – resistance out of range (the last digit is the zone number).
1520	A/D converter failed; incubator PCB defective.
1600	Computer control assay definition error. This error will occur for the following definitions: Well set, Wave scan, Checksum at the protocol sent from computer, Plate geometry, Filter(s), Features available, Mono
1700	Kinetic interval too short for selected options, or kinetic interval = 0. This error indicates that the kinetic interval in the current assay is too short. Increase the kinetic interval.
1900	Memory allocation failed. This error is typically used only for software development purposes. If it occurs, however, try turning the instrument off and then on again after a wait of 30 seconds. If the error persists, contact BioTek Technical Support.
1C00	A/D calibration standby signal on the analog board never went low, or A/D calibration standby line went low but never transitioned to a high. This error indicates there is a failure with the absorbance analog board initialization, or the cable to the PCB is defective.
2000	Barcode scanner did not see 10 characters from barcode. • Barcode positioned incorrectly on plate. • Insufficient carbon black in barcode label. • Barcode label not in Code 39 format.
2100	Invalid parameter value selected. This error can occur only during computer control, indicating that one of the following invalid assay configurations was sent to the instrument: • Temperature out of range • Wavelength not in ASCII format • Incorrect plate geometry • Incorrect row range or column range selected • Kinetic interval = 0 • Incorrect range or order selected for start wavelength / end wavelength • Incorrect well selected for scanning
2400	Middle sensor position incorrect.
	This error occurs when homing to the middle sensor and the optical flag is in a different position since the last autocalibration was performed.
2500 2502	Sweep mode read missed well location; last digit is the motor number.
2800- 2803	Motor currently in use; last digit indicates motor.

Code	Description and Probable Causes
2F00	Results data being is sent not acknowledged by host PC.
	This error indicates the handshaking between the host PC software and the reader did not complete. This is a lost data condition.
3200- 3201	Never saw A/D ready transition.
	This error indicates there is a failure with the absorbance analog board or the cable to the PCB is defective.

Fatal Errors

Fatal errors indicate conditions that require immediate attention. If a fatal error is displayed, contact BioTek's Technical Assistance Center for further instructions.

Code	Description
A100	Task control block not available.
A200	Read already in progress.
A300	Motors not available.
A301	Real time clock not available.
A302	Display not available.
A303	Flash not available.
A400	Failed code checksum test on power-up.
A600	Data flash write timed out.
A700	Data flash readback did not match write.
A800	Code flash write timed out.
A900	Memory allocation heap corrupted.

Instrument Dimensions

This chapter contains the PowerWave's dimensions, for use with robotic interfaces.

Instrument Dimensions

The figure below shows the location of the microplate carrier in reference to the exterior surfaces of the PowerWave, and the mounting holes on the bottom of the reader. This should facilitate system setup with a robotic unit.

If you purchased the BioStack to operate with the PowerWave, alignment hardware is included for positioning the instruments. For more information, refer to the BioStack Operator's Manual.

