Surface Areas of Cell Culture Vessels

The following are the surface areas of various tissue culture dishes or plates and their relative size to a 35 mm dish:

Dish Diameter	Actual Area	Calculated Area	Relative to 35 mm Dish
35 mm	10 cm ²	9.6 cm ²	1
60 mm	20 cm ²	28 cm ²	2
100 mm	60 cm ²	78.5 cm ²	6

Actual Area: These measurements were determined from actual measurements of common tissue culture vessels and the areas were rounded off.

Calculated Area: These values were calculated using the name of the culture vessel as the diameter.

Relative to 35 mm dish: This was calculated based on the actual area.

Area of a circle: The formula for calculating the area of a circle is: Area = $(pi)(r^2)$ where pi = 3.1416

:

Growth areas and yields of cultured cells* in various culture vessels			
Cell culture vessel	Growth area** (actual) in cm ²	No. of cells***	
Multiwell plates			
96	0.3	4-5 X 10 ⁴	
48	0.7	1.3 X 10 ⁵	
24	2	2.5 X 10 ⁵	
12	4	5.0 X 10 ⁵	
6	10	1.2 X 10 ⁶	
Dishes	-		
35 mm	10	1.2 X 10 ⁶	
60 mm	20	2.5 X 10 ⁶	
100 mm	60	7.5 X 10 ⁶	
145-150 mm	140	1.8 x 10 ⁷	
Flasks			
T25	25	3.1 X 10 ⁶	
T75	75	9.4 X 10 ⁶	
T150	150	1.9 X 10 ⁷	
T162	162	2.0 X 10 ⁷	
T165	165	2.1 X 10 ⁷	
40-50 ml	25	3.1 X 10 ⁶	
250-300 ml	75	9.4 X 10 ⁶	
650-750 ml	162-175	2.0 X 10 ⁷	
900 ml	225	2.8 X 10 ⁷	

*Such as NIH3T3, HeLa, CHO cells.

**Per well, can vary depending on the supplier. These values were determined from actual measurements of common tissue culture vessels and the values were rounded off.

***Assuming confluent growth. (~ $1.25 \times 10^5 \text{ cells/cm}^2$).