

TRU-SWEEP

Model 2800HT Ultrasonic Cleaner

Crest Ultrasonics Corp, USA

This Crest Tru-Sweep 2800 Heated Ultrasonic Cleaner has a tank capacity of 29½L (7.75Gal). This is a heated model.

Includes inner tray/basket to keep items from touching the bottom & sides. Lid is included.

Dimensions: 18"L x 12"W x 8"deep (inside) and 21"W x 13"D x 16"H (outside).

Inside Tank Dimensions: 19.5" x 11.5" x 8"

Overall Unit Dimensions: 20.3" x 12.15" x 15"

Tank Capacity: 7.75 Gallons (29.75 liters)

Internal thermostatically controlled heater.

All Brushed Stainless steel enclosures resist damage from chemical spills and will not crack.

Analog model with Zero to 30 minute mechanical timer.

Timer has continuous-operation capabilities.