

Blue 470nm LED light Transilluminator table

Specifications

BluView Transilluminator 470nm blue LED light

Window imaging size: 153x153mm

Synthetic paper hood with 580nm Amber filter
to picture via Smart phone

Bonus: Amber goggles (580nm) to cut bands.

Model: **CSMBE-200A**

Manufacture **Major Science, USA**



מפרט טכני:

Cat. No.	MBE-200A
Dimension (WxLxH)	Approx. 7.9"x7.9"x0.6" (200x200x13.9mm)
Viewing Area (WxL)	6"x6" (153x153mm)
Blue Light Source	15W
Illuminator Base Design	Flat Bed
Blue Light Wavelength	470nm
Automatic Shutdown	Approx. 6 min
Material	Aluminum Alloy
Power	DC 12V, 2A
Weight	Approx. 1.68lb (760g)