



Key Features

- Drawer Structure
- Strong Refrigeration
- Two Alarm Functions

More Detail:

http://www.haiermedical.com/products/preservation/40_deep_freezer/dw-40l92.html

Specifications

Technical Data

Type	-40 Deep Freeer
Cabinet Type	Upright
Climate Class	N
Cooling Type	Direct cooling
Defrost Mode	Manual
Refrigerant	CFC-Free
Noise Lever (dB)	40

Performance

Cooling performance()	-40
Temperature Range()	-10~-40

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50
--------------------	------------

Power(W)	220
Electrical Current(A)	2.1
Dimensions	
Capacity(L/Cu.Ft)	92/3.3
Net/Gross Weight(approx)	46/51(kg) 101.4/112.4(lbs)
Interior Dimensions(W*D*H)	465*440*645(mm) 18.3*17.3*25.4(in)
Exterior Dimensions(W*D*H)	597*610*860(mm) 23.5*24.0*33.9(in)
Packing Dimensions(W*D*H)	659*699*934(mm) 25.9*27.5*36.8(in)
Container load (20'/40'/40'H)	54/108/108
Functions	
High/Low Temperature	Y
Sensor Error	Y
Accessories	
Caster	Y
Foot	Y
Test Hole	N/A
Drawers	3
Temperature Recorder	N/A
USB Interface	N/A
Others	
Certificate	CE