


Decapitator



Stoelting's small animal decapitator has hardened and ground surgical stainless steel blades, aluminum base and stainless steel hardware. With its 1 1/2" square maximum blade opening, it is suitable for decapitating rats, mice and other small animals. The decapitator cuts cleanly and quickly through bone and tissue. The base has four holes for bench mounting.

The large animal decapitator is similar to our small animal decapitator, but with maximum blade opening of 4 inches square. For use with rabbits, small monkeys and other large sized animals.

Ordering Information:

51330 Small Animal Decapitator

51331 Large Animal Decapitator

51327 Replacement Blades for 51330 (FACTORY INSTALLED)