

מפרט טכני:

1. Fully stainless steel structure.
2. Indicator light indicates working state.
3. Over temperature&over pressure auto-protection.
4. Double scale indication pressure gauge.
5. Safe protection of water lacking.
6. Automatically shut off with beep reminding after sterilization.
7. Easy to operate,safe and reliable.
8. With two stainless steel sterilizing baskets.

▲ Specification:

Model	LS-B75L
Valid volume	75L(φ400*600 mm)
working pressure	0.22MPa
Working temperature	134 Deg.C
Max working pressure	0.23 Mpa
Heat average	$\leq \pm 1$ Deg.C
Timer	0~60min
Adjustment of temperature	105~134 Deg.C
Power	4.5KW /AC220V 50Hz
Overall dimension	560*560*980(mm)
Transport dimension	650*630*1150(mm)
N/G. W:	90Kg/70Kg

