

מדיח מעבדתי עם חימום דגם SMEG Laboratory dish washer GW4060S

GW4060S - LABORATORY GLASSWARE WASHER WITH DRYING SYSTEM

- washing up to 95°C - antibacterial thermal disinfection
- programs: 9 default + 6 free
- microprocessor programmable electronic control
- serial port RS232 for connection to external printer or PC
- forced hot air drying system with air filter 98% DOP
- peristaltic pump for liquid detergent
- peristaltic pump for acid liquid neutralizer
- possibility of installing one additional peristaltic pump
- washing chamber and internal door in AISI 316L stainless steel
- external panels made in AISI 304 stainless steel
- 2 independent washing levels through telescopic guides
- external dimensions: LxDxH: 600 x 600 x 850 mm
- internal working dimensions LxDxH: 525 x 490 x 570 mm
- incorporated water softener with salt sensor
- cold water and demineralized water connections
- electronic door lock system for user safety
- electrical connection: 3/N/PE 400 V - 50 Hz - 7 kW max
- single phase version available: 1/N/PE 230 V - 50 Hz - 2,8 kW max (GW4060S1)
- without trolleys



Machine optionals

- AF4060 - ABSOLUTE HEPA FILTER
- AS4060 - AQUA STOP WATER LEAK KIT
- B6040L - STAINLESS STEEL PLINTH FOR 60 CM MACHINES WITH LOCK
- FLUX4060 - DETERGENT DOSAGE CONTROL FOR ADDITIONAL PERISTALTIC PUMPS
- IRCP - CEILING SPRAY ARM
- P14060 - PERISTALTIC PUMP FOR LIQUID DETERGENT
- P34060 - PERISTALTIC PUMP FOR LIQUID ADDITIVES
- PAD1 - BOOSTER PUMP FOR NON-PRESSURE DEMINERALIZED WATER
- PAD2 - BOOSTER PUMP FOR NON-PRESSURE DEMINERALIZED WATER
- PAD2R - PAD2 RELAY KIT
- T6040 - STAINLESS STEEL FRAME FOR 60 CM MACHINES
- TOP60I - BUILT-IN TOP 60 CM
- WD-LANE - DATA COMMUNICATION CARD
- WD-LS3060 - DETERGENT LEVEL SENSOR
- WD-PRINTE - EXTERNAL PRINTER

סוגי מגרות אופציונליות בתשלום למדיח מעבדתי: SMEG Laboratory dish washer:

LM40SDS

UPPER FLASK WASHING WITH DRYING SYSTEM - 40 POSITIONS

- suitable for washing and drying of narrow neck glassware up to 225mm in height
- equipped with 40 nozzles U4140
- positioned on the upper level
- made of stainless steel



CS1-1

UPPER TROLLEY WITH SPRAYER

- suitable for holding of supports and baskets
- adjustable in two different positions (heights)
- made of stainless steel



CS2

LOWER TROLLEY

- suitable for holding of supports and baskets
- made of stainless steel



D-PF2

COMPLETE STAINLESS STEEL BOTTOM FOR LOWER TROLLEY

- made of stainless steel netting - mesh = 10 x 10 mm
- dim.: l x w x h = 305 x 425 x 11 mm
- positioned on lower trolley D-CS2



CP192

TEST TUBE BASKET - H = 165 MM

- capacity 160 standard test tubes
- equipped with 2 removable compartments with lid
- occupies 1/2 trolley D-CS2
- occupies 1/4 trolley CS1-1 / cs2
- made of stainless steel



CSK2

BASKET WITH HANDLES

- suitable for holding medium-sized instruments
- dimensions l x w x h 450 x 225 x 50 mm
- 5 x 5 mm mesh
- made of stainless steel



סוגי מגרות אופציונליות בתשלום למדיח מעבדתי: SMEG Laboratory dish washer (המשך)

LPM2020DS - MIXED WASHING TROLLEY FLASKS / PIPETTES – 20 + 20 POSITIONS WITH DRYING SYSTEM

- suitable for washing and drying 20 pipettes up to 55 cm + 20 pieces of narrow-neck glassware equipped with the following nozzles :
- equipped with the following nozzles : 2 x U6260, 2 x U6240, 2 x U6220, 3 x U4180, 2 x U4160, 2 x U4140, 3 x U3110, 4 x U390, 1 x UC6, 20 x ULB40 for pipettes max height 55 cm
- positioned on lower level
- made of stainless steel



LM80DS - UNIVERSAL FLASK WASHING WITH DRYING SYSTEM - TWO WASHING LEVELS

- suitable for washing and drying of narrow neck glassware up to 225mm in height
- capacity 68 pieces
- equipped with 68 nozzles U4140 + 1 nozzle UC6
- positioned on 2 independent levels
- made of stainless steel



B6040L - STAINLESS STEEL PLINTH FOR 60 CM MACHINES WITH LOCK

- dimensions: l x w x h = 600 x 530 x 400 mm
- rises the machine working level to 70 cm from the floor
- detergent cabinet with front door
- equipped with adjustable feet
- detergent cabinet with anti-drip tray
- front door with lock
- made in stainless steel

